## L-NEW 8 STEPS

## For Allard AFO Customization

## "Get The Alignment Right"


Product Selection. Decide stability, based on patient needs



- **Size Selection Correct Application, based** on the following:
  - Foot length
- Activity level
- Body type
- Spasticity
- Deficit type
- Foot Orthotic (custom or prefab), UCBL, or SMO



**Shoe Selection & Heel Height** 



 Accommodate heel height difference!

#### **Frontal Plane Alignment**

- With the patient standing safely in parallel bars, check their alignment with either a plumb line or square.
- Make accommodation on plantar surface of the product using crepe or cork.



#### **Proximal Control**

Add as little as 1/16" (1.6mm) heel wedge



#### **Patient Comfort**

- Slide foot plate forward or back in the shoe to distribute the pressure evenly across the tibia.
- Add SoftKIT or ComfortKIT.



### **Patient Education.**

#### DONT's

- Use hand support ON Sit-to-Stand
- Stairs Ball of Foot
- SQUATS not at all
- Use SoftKIT, ComfortKIT or Custom Interface for tibia as provided with your brace.
- · Use interface on footplate
- Inspect brace on a daily basis
- Check SKIN on a daily basis
- Wear recommended footwear

HEEL HEIGHTS

KiddieGAIT® Baby KiddieGAIT® S/M=1/4"(7mm)

L=5/16"(8mm)

KiddieROCKER® S=5/16"(8mm) M/L/XL=3/8"(9mm) L/XL=5/8"(16mm)

ToeOFF®/BlueROCKER® Ypsilon® FLOW 1/2 XS/S=7/16"(12mm) M=9/16"(15mm)

All 2 1/2 & SHORT models XS-XL=7mm



# Step by Step Guide for ToeOFF®/BlueROCKER®

## **Fitting & Product Selection**

Evaluate Patient with standard ToeOFF®/BlueROCKER®.

Establishes a baseline for altering the dynamic properties of a standard product.

Choose ToeOFF® or BlueROCKER® based on Patient's height, weight, shoe size and proximal deficits.

Choose proper shoe.

Measure Heel Height of Shoe.

Accommodate ToeOFF®/BlueROCKER®
with heel/toe lift if necessary to match
heel rise of ToeOFF®/BlueROCKER®.



#### **Shoe Selection**

- Shank & counter
- Toe-to-heel
- Rocker sole
- Facilitate 3rd rocker

Ypalor\* TooOff\* BluRDCKER\*\*

2 Capture the Patient model with Cast, Scan or Measurement

#### Cast:

- I. Posterior/Medial Opening
- 2. Anterior & Lateral Plumb lines

#### Scan:


File extension type: \*.obj , .ply, .stl, .vsrf, .cxpxp, .cln, .ssn and Vorum Spectra scan.

#### **Measurements:**

- 1. Fill out pages5-7 Completely
- 2. Picture required.

Fill out Custom Forms Completely

Measure Heel Height of Shoe: Critical to get proper alignment of ToeOFF or BlueROCKER. Either measure with a ruler or with calipers.







On order form - Address "Change M-L Stiffness", "Change A-P Stiffness", "Relief of Tibial Crest" and "Relief of OTHER Bony Prominences"

Easy as 1,2,3

Make sure custom forms are filled out properly!

