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**8.** Firm heel counter – for proper control of the rear foot  
 • The foot and the orthosis should be secured by the shoe. Be aware of high tone.  
 • Rocker-type sole at toe end for smoother toe-off.  
 • Laced - for easier donning and doffing, and to allow adjustable compressive support over the mid-foot.  
 • Rubber sole – to minimise the chance of slipping on wet surfaces.

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**12.** • Removable insole – to allow space for the footplate with modifications.  
 • Position in the shoe. The strut should be located just posterior to the 5th metatarsal head and extend upward without touching the tibia. Shift the footplate forward or backward to achieve this proper alignment. Trim or add material to the footplate to prevent it from sliding in the shoe and losing proper alignment.  
 a) Add material medially to keep 5th MTP free from pressure.  
 b) Add material laterally if too far from 5th MTP.  
 c) Shift the footplate forward or backwards to correct knee position and to avoid contact with tibia crest.

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**GB** Thermoset Hybrid Compound  
 • Do not heat.  
 • Do not drill into material

**11.** Use the sizing chart as a guide only. As the size increases, the rigidity and the height of the orthosis increases. Consider sizing up to increase strength of the orthosis or down to increase flexibility. Refer to ToeOFF™ Clinical Manual for more information.

**12.** • Firm heel counter – for proper control of the rear foot  
 • The foot and the orthosis should be secured by the shoe. Be aware of high tone.  
 • Rocker-type sole at toe end for smoother toe-off.  
 • Laced - for easier donning and doffing, and to allow adjustable compressive support over the mid-foot.  
 • Rubber sole – to minimise the chance of slipping on wet surfaces.

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GB	<p>All three products (ToeOFF®, Ypsilon™ and BlueRocker™) are fabricated to accommodate approximately 1.5 cm (5/8 inch) heel height. If the shoe has a significantly higher or lower heel, add or reduce the effective heel height. To encourage: Knee Extension (minimize flexion) → Decrease Heel Height. Knee Flexion (delay extension) → Increase Heel Height. To influence towards knee valgus → wedge the bottom of the footplate on the lateral side. To influence towards knee varus → wedge the bottom of the footplate on the medial side.</p>	<p>Pressure at Calfband - Ypsilon™ a) Angle the strap to fit the shape of the calf. b) Add 3 mm soft foam padding to the strap.</p>	<p>Patient Comfort A. The orthosis should always have padding on the inside before delivered to the patient. If pressure on the tibial crest add or reduce both laterally and medially, leaving an open channel for tibial crest relief. B. Pressure Over Instep. Shape and adhere padding to the distal medial aspect of the lateral strut to shift the foot medially. See instructions # 5. C. Pressure at Calfband. Add 3mm soft foam padding to strap. SoftKIT™ and SoftShell™ are two ready made options for easy, simple and quick padding.</p>	<p>To optimise gait and maximise product durability, correct the position of the foot to achieve a position as close to neutral as possible when bearing weight. The most common foot related problems like pronation, supination, pes varus, pes valgus should be corrected with an insole. Accommodation of other foot orthotics: use rubber cement to adhere metatarsal pads, arch supports, pads for posting etc. to the top of the footplate. To increase patient comfort - always cover the footplate with leather or similar material.</p>	<p>Be aware that the orthosis may alter the patient's abilities for some activities. I.e., driving a car. Always make sure the patient reads the patient instructions supplied with the products. Inform the patient that some pain can occur in muscles that have not been active when starting to use a new orthotic.</p>	<p>More detailed information about fitting is available on the ToeOFF®/Ypsilon™ CD-Rom, or at www.campscandinavia.se. Please read the ToeOFF® Gait Management Manual for complete instructions.</p>		
----	--	---	---	---	---	---	--	--

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YPSILON® ToeOFF® BLUEROCKER™

Manufactured by:  
Camp Scandinavia AB, Sweden  
Phone: +46 42 25 27 01  
www.campscandinavia.se

ToeOFF®, Ypsilon™ and BlueRocker™ are patented products.  
© Camp Scandinavia AB / aug 05

- Utprovninginstruktion
- Tilpassingsmanual
- Tilpassingsmanual
- Sovitusohjeet
- Professional Instructions
- Technische Anleitung
- Professionele instructies
- Manuale Professionale
- Manual Profesional
- Manual Clinico
- Manuel Professionnel