YPSILON™

Utprovninginstruktion

Tilpasningsmanual

Utprøvningsmanual

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Professional Instruction

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Professionele instructies

Manual Professional

Manual Profesional

Manual Clinico

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2.

When Grinding/Cutting

- Protect the eyes
- Cover the nose and mouth.
- Do not grind/cut into Carbon or Kevlar
- Do not overheat the composite

3.

Patient Assessment

Before choosing product style or size, make a proper assessment of the affected foot to determine function and stability in both open and closed chain. Take into consideration the larger the size, the more rigidity. When assessment shows ankle instability, ToeOFF® is recommended instead of YPSILON™. More about patient assessment is described in the ToeOFF® Clinical Manual.

4.

- Firm heel counter – for proper control of the rear foot
- Rocker-type sole at toe end for smoother toe-off
- Laced - for easier donning and doffing, and to allow adjustable compressive support over the mid-foot.
- Rubber sole – to minimise the chance of slipping on wet surfaces.
- Removable insole – to allow space for the YPSILON™ footplate with modifications.
- 1,5 cm (5/8 inches) toe-to-heel height differential – as a starting point to control the knee extension moment.

5.

Use the sizing chart as a guide only. As the size of the YPSILON™ increases, the *rigidity* and the *height* of the orthosis increase. Consider sizing up to increase strength of the orthosis or down to increase flexibility. **Refer to ToeOFF® Clinical Manual for more information.**

6.

Instep clearance

Position YPSILON™ in the shoe. The strut should be located just posterior to the 5th metatarsal head and extend upward without touching the tibia. Shift the footplate forward or backward to achieve this proper alignment. Trim or add material to the footplate to prevent it from sliding in the shoe and losing proper alignment.

- a) Shift the footplate lateral to keep 5th MTP free from pressure.
- b) Shift the footplate medially if too far from 5th MTP.
- c) Shift the footplate forward or backwards to correct position and to avoid contact with tibia crest.

7.

Pressure at Calfband

- a) Angle the strap to fit the shape of the calf.
- b) Add 3mm soft foam padding to the strap.

8.

To optimise gait and maximise product durability, correct the position of the foot to achieve a position as close to neutral as possible when bearing weight.

Pronation – Providing there are no other foot deformities, post the medial aspect of the calcaneus on the footplate to decelerate the pronation moment and vice versa for supination.

Accommodation of other foot orthotics. Use rubber cement to adhere metatarsal pads, arch supports, pads for posting etc. to the top of the footplate.

9.

YPSILON™ is fabricated to accommodate approximately 1,5 cm (5/8 inch) heel height. If the shoe has a significantly higher or lower heel, add or reduce the heel height.

To encourage:

Knee Extension (minimize flexion) → Decrease Heel Height

Knee Flexion (delay extension) → Increase Heel Height

Knee Valgus → Wedge bottom of footplate on lateral side

Knee Varus → Wedge bottom of footplate on medial side

10.

Please read the "ToeOFF® Clinical Manual" for complete instructions.

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