
allard | ^N_T

Carbon Footplate

28035

Kolfibersula / Kulfibersål / Kullfibersåle
Hiilikuitupohjainen

Instructions for use
Användarinstruktion
Brugsanvisning
Anvendelseinstruktion
Käyttöohjeet
Gebrauchsanweisung
Gebruiksaanwijzing
Instructions d'utilisation
Istruzioni per l'uso
Instrucciones de uso

ENGLISH

Product description

The carbon footplate limits the range-of motion (ROM) in the metatarsophalangeal joints MTP-joints, distributes the pressure in the forefoot and facilitates roll-over.

Intended use

The footplate is intended to protect the metatarsophalangeal joints MTP-joints and the forefoot by limiting the range of motion (ROM), distribute the pressure on the forefoot and facilitate roll-over. The footplate can also be used after immobilizing treatment with a POST-OP shoe or walker in order to gradually increase movement, and to improve balance after toe amputation of digits IV and/or V. Must be covered by an insole. Single patient, multiple use.

Indications

Hallux rigidus, hallux limitus, stress fractures, forefoot trauma, toe sprain and other injuries in the forefoot/midfoot.

Contraindications

Open ulcers on the foot, non stable fractures and a weight over 300 lbs/ 135 kg.

Adjustments

This product should only be altered by a medical professional.

Warnings

Extra precaution is needed if the patient has neuropathy, diabetes or other conditions that may lead to ulcers. Make sure the carbon footplate is covered by a functional foot orthotic that manages the plantar surface and the selected shoe has sufficient inner volume so pressure sores do not occur.

Fiber reinforced thermoset composite: Do not heat the material

When grinding/cutting: Protect the eyes, use dust protection over nose and mouth, make sure that the edges of the footplate are carefully smoothed and all loose fiber particles are removed. Do not overheat (max. 100° C/212° F) when grinding the composite.

Material Specifications

Carbon fiber in an epoxy matrix.

User instruction

The condition of the foot should be checked daily. In the case of skin irritation or discoloration, discontinue use immediately and contact the fitter/prescriber. If your ability to monitor the foot daily is reduced, ask for assistance. Do not bend the carbon footplate excessively, for example, by squatting.

Washing Instructions

Wash with soap and water.


Disposal

To be disposed in accordance with local policies and regulations.

Information


Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in the location where the user and/or patient is established.

Adjusting the carbon footplate

Remove the insole from the shoe.


Choose the carbon footplate closest to the patient's shoe size.


Use the insole from the footwear as a template. Trace the template onto the carbon footplate.

The insole should be a couple of millimeters larger than the carbon footplate.

Grind the carbon footplate with a belt sander or Dremel and smooth the edges.

Place the carbon footplate in the shoe.

Return the shoe insole on top of the carbon footplate.

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Camp Scandinavia AB, Karbingatan 38
Camp Scandinavia AB