

PARTIAL FOOT Order Form

All Partial Foot orders must be shipped to our Chattanooga location.

Date: _____

Contact: _____	SHIP VIA:	<input type="radio"/> UPS	<input type="radio"/> FedEx	<input type="radio"/> Cust. P/U
P.O. #: _____				
Phone: _____	SERVICE:	<input type="radio"/> Ground	<input type="radio"/> 3-Day	<input type="radio"/> 2-Day
Fax: _____		<input type="radio"/> ND Saver	<input type="radio"/> Next Day	<input type="radio"/> ND Early AM
Email: _____				

SHIP TO:	BILL TO:	<input type="radio"/> Same as Ship To
Acct #: _____	Acct #: _____	
Name: _____	Name: _____	
Address: _____	Address: _____	
Address 2: _____	Address 2: _____	
City, State, Zip Code: _____	City, State, Zip Code: _____	

PATIENT INFO

Patient ID: _____

Height: _____ Weight: _____

SHOE INFO

Included shoes with cast (Recommended) Shoe insert tracing Shoe Size: _____ Sound Side Foot Length: _____

AMPUTATION SIDE

Left Right

AMPUTATION LEVEL

Great Toe TMA Lisfranc Chopart

Other _____ Surgical Congenital

Date of Amputation: _____

ALLARD AFO CHOICE

BlueROCKER® 2½ Recommended for all PFAs proximal to and including great toe amputations

ToeOFF® 2½ May be considered for 2 - 5 toe amputations

PADDING OPTIONS

ComfortKIT™ (Recommended) Premium 5mm memory foam interface for patients requiring added cushioning and a more intimate fit due to the shape of the tibia, skin condition, diabetes, and/or activity level.

SoftKIT™ Soft foam in 3mm.

CoverKIT™ Soft interface combined with a spandex sleeve to cover the AFO (Recommended for venous insufficiency)

DISTAL CUSHION

Yes


No

Cont'd

CARBON FOOT PLATE ONLY (Requires tracing of foot or shoe insert)

- Left Qty _____
- Right Qty _____

STANDING LIMB LENGTH - Fibular head to floor



SAGITTAL PLANE CALCANEAL INCLINATION ANGLE

Weightbearing Actual

- 25°
- 20°
- 15°
- 10°
- 5°
- 0°



Available ROM

- 25°
- 20°
- 15°
- 10°
- 5°
- 0°




Wedge Fabrication Order

- 25°
- 20°
- 15°
- 10°
- 5°
- 0°



FRONTAL PLANE CALCANEAL INVERSION-EVERSION

Weightbearing Actual



+6 +4 +2 0 -2 -4 -6 +6 +4 +2 0 -2 -4 -6


circle one

circle one

Inversion Eversion

Eversion Inversion

Post-Fabrication Order



+6 +4 +2 0 -2 -4 -6 +6 +4 +2 0 -2 -4 -6

circle one

circle one

Lateral Medial

Medial Lateral

SPECIAL INSTRUCTIONS:
