


allard|USA



*comfort*KIT™

**A new soft and self-molding interface
option for ToeOFF® and BlueRocker™**

A TOEOFF® FAMILY PRODUCT

comfortKIT™

A new soft solution



Let's face it. Everyone's leg has a different shape and different skin condition. The new comfortKIT offers a solution for ToeOFF & BlueRocker wearers who due to the shape of the tibia, skin condition, and/or activity level, will benefit from an interface that offers added cushioning and more intimately contours to the shape of the leg than the standard SoftKIT.

An inner layer of polyurethane foam molds naturally to the leg, yet offers the "bounce" to return to its original shape when not in use.

Like SoftKIT, comfortKIT is easily adhered to ToeOFF or BlueRocker with MikroFix aggressive velcro kit.

The comfortKIT consists of two interfaces and one set of MikroFix velcro kit.

Item No.	Size
28398	XSmall, Small, Medium, Large/XLarge

