

CSUS


by Allard USA

C|s

SHOULDER STABILIZER
Rx Only — Single Patient Use Only

C|s

Instructions for Use



FEEL THE FUNCTION

The intended function of this stabilizer is for treatment of traumatic or post-surgical immobilization of the shoulder. Review instructions on back page carefully before application.

For questions, please contact your physician, the professional who fit this device or Allard USA Customer Service:

CALL: 888-678-6548 FAX: 800-289-0809
info@allardusa.com

Allard USA, Inc. 300 Forge Way, Suite 3 Rockaway, NJ 07866-2056
www.allardusa.com

CSUS

by Allard USA

C|s

SHOULDER STABILIZER
Rx Only — Single Patient Use Only

C|s

Instructions for Use



FEEL THE FUNCTION

The intended function of this stabilizer is for treatment of traumatic or post-surgical immobilization of the shoulder. Review instructions on back page carefully before application.



For questions, please contact your physician, the professional who fit this device or Allard USA Customer Service:

CALL: 888-678-6548 FAX: 800-289-0809
info@allardusa.com

Allard USA, Inc. 300 Forge Way, Suite 3 Rockaway, NJ 07866-2056
www.allardusa.com

Shoulder Stabilizer

1. Wrap the support around the rib cage, positioning between the sternum and the waist, with the attached arm band on the side of the injured arm. Wrap the attached arm band around the affected arm right above the elbow and secure.



2. Wrap the small band around the affected wrist and secure in clinician recommended location on front of large panel.



3. Pass the narrow strap over the unaffected shoulder, feed through both D-Rings and adjust length to desired position.

Indications / Intended Use

Shoulder injuries and post surgical immobilization.

Contraindications

Inability to bend arm to at least 90°.

Wash instruction

Hand wash in cold water using mild detergent.
Air dry.

Adjustment



This stabilizer should only be fitted by a licensed medical professional.

Warning

If wearing the stabilizer causes any unusual reactions, swelling or prolonged pain, remove and consult your medical professional. Check skin on a daily basis.

Shoulder Stabilizer

1. Wrap the support around the rib cage, positioning between the sternum and the waist, with the attached arm band on the side of the injured arm. Wrap the attached arm band around the affected arm right above the elbow and secure.



2. Wrap the small band around the affected wrist and secure in clinician recommended location on front of large panel.



3. Pass the narrow strap over the unaffected shoulder, feed through both D-Rings and adjust length to desired position.

Indications / Intended Use

Shoulder injuries and post surgical immobilization.

Contraindications

Inability to bend arm to at least 90°.

Wash instruction

Hand wash in cold water using mild detergent.
Air dry.

Adjustment

This stabilizer should only be fitted by a licensed medical professional.

Warning

If wearing the stabilizer causes any unusual reactions, swelling or prolonged pain, remove and consult your medical professional. Check skin on a daily basis.