

CSUS

by Allard USA

C|s

SLING & SWATHE
Rx Only — Single Patient Use Only

C|s

Instructions for Use



FEEL THE FUNCTION

The intended function of this sling is to transfer the weight of the forearm to the patient's back and shoulders with the option to secure the affected arm against the body. Review instructions on back page carefully before application.

For questions, please contact your physician, the professional who fit this device or Allard USA Customer Service:

CALL: 888-678-6548 FAX: 800-289-0809
info@allardusa.com

Allard USA, Inc. 300 Forge Way, Suite 3 Rockaway, NJ 07866-2056
www.allardusa.com

CSUS

by Allard USA

C|s

SLING & SWATHE
Rx Only — Single Patient Use Only

C|s

Instructions for Use



FEEL THE FUNCTION

The intended function of this sling is to transfer the weight of the forearm to the patient's back and shoulders with the option to secure the affected arm against the body. Review instructions on back page carefully before application.

For questions, please contact your physician, the professional who fit this device or Allard USA Customer Service:

CALL: 888-678-6548 FAX: 800-289-0809
info@allardusa.com

Allard USA, Inc. 300 Forge Way, Suite 3 Rockaway, NJ 07866-2056
www.allardusa.com

Sling & Swathe

1. Place affected arm inside the sling. Raise buckle up and slip straps over head.



2. Wrap strap ends under each arm and secure to front of sling. Adjust straps as needed.

3. Wrap the swathe over the affected arm and under the opposite arm. Secure for snug, comfortable fit.



Indications / Intended Use

Shoulder injuries and post surgical immobilization.

Contraindications

Inability to bend arm to at least 90°.

Wash instruction

Hand wash in cold water using mild detergent.
Air dry.

Adjustment



This sling should only be fitted by a licensed medical professional.

Warning

If wearing the sling causes any unusual reactions, swelling or prolonged pain, remove and consult your medical professional. Check skin on a daily basis.


Sling & Swathe

1. Place affected arm inside the sling. Raise buckle up and slip straps over head.



2. Wrap strap ends under each arm and secure to front of sling. Adjust straps as needed.

3. Wrap the swathe over the affected arm and under the opposite arm. Secure for snug, comfortable fit.



Indications / Intended Use

Shoulder injuries and post surgical immobilization.

Contraindications

Inability to bend arm to at least 90°.

Wash instruction

Hand wash in cold water using mild detergent.
Air dry.

Adjustment

This sling should only be fitted by a licensed medical professional.

Warning

If wearing the sling causes any unusual reactions, swelling or prolonged pain, remove and consult your medical professional. Check skin on a daily basis.