

allard|USA

CENTRAL FABRICATION  
SUPPLIES & ACCESSORIES

[allardusa.com](http://allardusa.com)






## CHAFE W/ D-RING

Chafes are made from 100% polyester material with a carbon steel, zinc plated d-ring and an aluminum eyelet. Choose Black or White. Available in 3 widths. Sold in packages of 25 pieces.

Item No.	Description	Size	Color
47489 0001	Chafe with 1" metal D-Ring	1"	Black
47489 0002	Chafe with 1.5" metal D-Ring	1.5"	Black
47489 0003	Chafe with 2" metal D-Ring	2"	Black
47491 0001	Chafe with 1" metal D-Ring	1"	White
47491 0002	Chafe with 1.5" metal D-Ring	1.5"	White
47491 0003	Chafe with 2" metal D-Ring	2"	White



## NYLON STRAP W/ EYELET

Heat welded hook and loop straps are 100% polyamide and include a reinforced aluminum eyelet. Choose Black or White. Available in 3 widths. Sold in packages of 25 pieces.



Item No.	Description	Size	Color
47492 0001	3" hook/9" loop	1" x 12"	Black
47492 0002	4" hook/12" loop	1" x 16"	Black
47492 0003	5" hook/15" loop	1" x 20"	Black
47492 0011	3" hook/9" loop	1.5" x 12"	Black
47492 0012	4" hook/12" loop	1.5" x 16"	Black
47492 0013	5" hook/15" loop	1.5" x 20"	Black
47492 0021	3" hook/9" loop	2" x 12"	Black
47492 0022	4" hook/12" loop	2" x 16"	Black
47492 0023	5" hook/15" loop	2" x 20"	Black
47494 0001	3" hook/9" loop	1" x 12"	White
47494 0002	4" hook/12" loop	1" x 16"	White
47494 0003	5" hook/15" loop	1" x 20"	White
47494 0011	3" hook/9" loop	1.5" x 12"	White
47494 0012	4" hook/12" loop	1.5" x 16"	White
47494 0013	5" hook/15" loop	1.5" x 20"	White
47494 0021	3" hook/9" loop	2" x 12"	White
47494 0022	4" hook/12" loop	2" x 16"	White
47494 0023	5" hook/15" loop	2" x 20"	White

## NYLON STRAP WITH CHAFE

Chafe made from 100% polyester, heat welded. Includes a carbon steel, zinc plated d-ring. Straps are made from 100% polyamide (nylon). Straps are heat welded hook and loop. Choose Black or White. Available in 3 widths. Sold in packages of 25 pieces.


Item No.	Description	Size	Color
47495 0001	3" hook/9" loop	1" x 12"	Black
47495 0002	4" hook/12" loop	1" x 16"	Black
47495 0003	5" hook/15" loop	1" x 20"	Black
47495 0011	3" hook/9" loop	1.5" x 12"	Black
47495 0012	4" hook/12" loop	1.5" x 16"	Black
47495 0013	5" hook/15" loop	1.5" x 20"	Black
47495 0021	3" hook/9" loop	2" x 12"	Black
47495 0022	4" hook/12" loop	2" x 16"	Black
47495 0023	5" hook/15" loop	2" x 20"	Black
47497 0001	3" hook/9" loop	1" x 12"	White
47497 0002	4" hook/12" loop	1" x 16"	White
47497 0003	5" hook/15" loop	1" x 20"	White
47497 0011	3" hook/9" loop	1.5" x 12"	White
47497 0012	4" hook/12" loop	1.5" x 16"	White
47497 0013	5" hook/15" loop	1.5" x 20"	White
47497 0021	3" hook/9" loop	2" x 12"	White
47497 0022	4" hook/12" loop	2" x 16"	White
47497 0023	5" hook/15" loop	2" x 20"	White





## DACRON® BACKED STRAP WITH EYELET

Hook and loop stitched to Dacron® material for strength. Includes a reinforced aluminum eyelet on one end. Score lines serve as guides to trim straps. Choose between 3 widths. Available in Black & White. Sold in packages of 25 pieces.

Item No.	Description	Size	Color
47485 0001	2" hook/4" loop	1" x 6"	Black
47485 0002	3" hook/9" loop	1" x 12"	Black
47485 0003	4" hook/12" loop	1" x 16"	Black
47485 0004	5" hook/13" loop	1.5" x 18"	Black
47485 0005	5" hook/15" loop	2" x 20"	Black
47487 0001	2" hook/4" loop	1" x 6"	White
47487 0002	3" hook/9" loop	1" x 12"	White
47487 0003	4" hook/12" loop	1" x 16"	White
47487 0004	5" hook/13" loop	1.5" x 18"	White
47487 0005	5" hook/15" loop	2" x 20"	White





## NYLON STRAP

Heat welded hook and loop straps are 100% polyamide (nylon). Choose Black or White. Available in 3 widths. Sold in packages of 25 pieces.

Item No.	Description	Size	Color
47488 0001	4" hook/8" loop	1" x 12"	Black
47488 0002	4" hook/12" loop	1" x 16"	Black
47488 0003	5" hook/15" loop	1" x 20"	Black
47488 0011	4" hook/8" loop	1.5" x 12"	Black
47488 0012	4" hook/12" loop	1.5" x 16"	Black
47488 0013	5" hook/15" loop	1.5" x 20"	Black
47488 0021	4" hook/8" loop	2" x 12"	Black
47488 0022	4" hook/12" loop	2" x 16"	Black
47488 0023	5" hook/15" loop	2" x 20"	Black
47486 0001	4" hook/8" loop	1" x 12"	White
47486 0002	4" hook/12" loop	1" x 16"	White
47486 0003	5" hook/15" loop	1" x 20"	White
47486 0011	4" hook/8" loop	1.5" x 12"	White
47486 0012	4" hook/12" loop	1.5" x 16"	White
47486 0013	5" hook/15" loop	1.5" x 20"	White
47486 0021	4" hook/8" loop	2" x 12"	White
47486 0022	4" hook/12" loop	2" x 16"	White
47486 0023	5" hook/15" loop	2" x 20"	White



## LEATHER STRAP W/ EYELET

Leather straps have an aluminum eyelet on one end and 100% polyester loop sewn on opposite side. Score lines serve as guides to trim straps without fraying. Only available in Dark Brown. Choose from 3 widths . Sold in packages of 25 pieces.

Item No.	Description	Size	Color
47499 0001	Leather strap w/ eyelet	1" x 4.5"	Dark Brown
47499 0002	Leather strap w/ eyelet	1" x 5.5"	Dark Brown
47499 0003	Leather strap w/ eyelet	1.5" x 6.5"	Dark Brown
47499 0004	Leather strap w/ eyelet	1.5" x 7.75"	Dark Brown
47499 0005	Leather strap w/ eyelet	2" x 11"	Dark Brown

## AFO STRAP PAD

AFO Strap Pads provide cushioning on a pressure sensitive area. Nylon and polyester plush loop fabric has 100% SBR incorporated into the layers for added relief on pressure points. Only available in Black. Choose from 2 widths. Sold in packages of 25 pieces.

Item No.	Description	Size	Color
47498 0001	AFO Strap Pad, fits 1/2" wide strap	1" x 1"	Black
47498 0002	AFO Strap Pad, fits 3/4" wide strap	1" x 1.5"	Black
47498 0003	AFO Strap Pad, fits 1 1/4" wide strap	1" x 2"	Black
47498 0021	AFO Strap Pad, fits 1/2" wide strap	2" x 1"	Black
47498 0022	AFO Strap Pad, fits 3/4" wide strap	2" x 1.5"	Black
47498 0023	AFO Strap Pad, fits 1 1/4" wide strap	2" x 2"	Black



## KLICK BUCKLES

Klick buckles are made from 100% polyamide material. For many applications but intended to go with appropriately sized Nylon Chafes, Alligator Tabs & Double Faced Loop. Only available in Black. Choose from 2 widths. Sold in packages of 25 sets.

Item No.	Description	Size	Color
69133 0001	Klick Buckles	1"	Black
69133 0002	Klick Buckles	1.5"	Black



## NYLON CHAFE

Nylon chafes are made from 100% polyamide material. Made for many applications but intended to go with appropriately sized Klick Buckles, Alligator Tabs & Double Faced Loop. Available in Black only. Choose from 2 widths. Sold in packages of 25 pieces.

Item No.	Description	Size	Color
69134 0001	Nylon Chafe	1"	Black
69134 0002	Nylon Chafe	1.5"	Black





## PADDED STRAP

Made from polyurethane & polyamide, this 1 ½” wide roll of padding is soft and plush. Without non-slip silicone on the back side. Only available in Black. Sold in packages of one 25 yard roll.

Item No.	Description	Size	Color
69135 0000	Padded Strap Roll	40 mm x 25 yds	Black



## ALLIGATOR TABS

Polyester alligator tabs are easy to use on our Double Faced Loop to create an instant strap without the need for sewing. Excellent option for patients with limited hand dexterity. Fits appropriately sized Klick Buckles, Chafes & Double Faced Loop. Only available in Black. Choose from 2 widths. Sold in packages of 10 pieces.

Item No.	Description	Size	Color
69136 0000	Alligator Tabs	16 mm wide	Black
69136 0002	Alligator Tabs	25 mm wide	Black



## VELCRO® DOUBLE FACED LOOP

This polyester Double Faced Loop is heat joined and will not separate. It is compatible with all conventional hook materials. Fits appropriately sized Klick Buckles, Chafes & Alligator Tabs. Choose from 2 widths. Only available in Black. Sold in packages of one 27 yard roll.

Item No.	Description	Size	Color
69137 0000	Velcro® Double Faced Loop Roll	25mm wide x 27 yds	Black
69137 0002	Velcro® Double Faced Loop Roll	38mm wide x 27 yds	Black



## ALLARD AFO STRAP SETS

Allard AFO straps are developed to meet the needs to fit from KiddieGAIT® Baby sizes to BlueROCKER®. Straps are available with or without D-Rings. Depending on hand dexterity, straps can be mounted on either side of the AFO. Only available in Black. Baby sizes sold in packages of one strap. Adult sizes sold in packages of 2 (upper and lower strap).

Item No.	Description	Size	Color
28373 0008	KiddieGAIT® Strap	Baby Small	Black
28373 0009	KiddieGAIT® Strap	Baby Medium	Black
28373 0010	KiddieGAIT® Strap	Baby Large	Black
28368 0011	KiddieGAIT® & KiddieROCKER®, No D-Ring	Small	Black
28368 0012	KiddieGAIT® & KiddieROCKER®, No D-Ring	Medium	Black
28368 0013	KiddieGAIT® & KiddieROCKER®, No D-Ring	Large/X-Large	Black
28369 0011	KiddieGAIT® & KiddieROCKER®, W/D-Ring	Small	Black
28369 0012	KiddieGAIT® & KiddieROCKER®, W/D-Ring	Medium	Black
28369 0013	KiddieGAIT® & KiddieROCKER®, W/D-Ring	Large/X-Large	Black
28760 0010	2.0 & 2½ Allard AFO Strap, No D-Ring	X-Small	Black
28760 0011	2.0 & 2½ Allard AFO Strap, No D-Ring	Small	Black
28760 0012	2.0 & 2½ Allard AFO Strap, No D-Ring	Medium	Black
28760 0013	2.0 & 2½ Allard AFO Strap, No D-Ring	Large	Black
28760 0014	2.0 & 2½ Allard AFO Strap, No D-Ring	X-Large	Black
28761 0010	2.0 & 2½ Allard AFO Strap, W/D-Ring	X-Small	Black
28761 0011	2.0 & 2½ Allard AFO Strap, W/D-Ring	Small	Black
28761 0012	2.0 & 2½ Allard AFO Strap, W/D-Ring	Medium	Black
28761 0013	2.0 & 2½ Allard AFO Strap, W/D-Ring	Large	Black
28761 0014	2.0 & 2½ Allard AFO Strap, W/D-Ring	X-Large	Black





## MIKROFIX™ FOR ALLARD AFO WINGS

MikroFIX™ Kits are precut shapes used to adhere straps to the Allard AFO. Self-adhesive backing for easy application. Use to replace MikroFIX™ on wings of Allard AFOs. MikroFIX™ Kit for 2.0 & 2 ½ Allard AFOs have precut shapes for wings only. MikroFIX™ Kit for Original Allard AFOs includes Tear Drop Shape for AFO Padding.

Item No.	Description	Size	Color
28396 0011	Original Allard AFOs	XS/S	Black
28396 0013	Original Allard AFOs	M/L/XL	Black
28754 0000	ToeOFF® & BlueROCKER® 2.0 & 2½	One Size	Black



## MIKROFIX™ FOR ALLARD AFO PADDING

MikroFIX™ Tear Drop Shape Set for 2.0 & 2 ½ AFOs. MikroFIX™ Kits are precut Velcro® shapes used to adhere the SoftKIT™ or ComfortKIT™ to the Allard AFO. Self-adhesive backing for easy application. Intended as replacements for the originals supplied with SoftKIT™ and ComfortKIT™. Only available in Black. Choose appropriate size for AFO being fit. Pediatric sizes sold in packages of 1 piece. Adult sizes sold in packages of 2 pieces.

Item No.	Description	Size	Color
28379 0011	MikroFIX™ for KiddieGAIT®/KiddieROCKER®	S	Black
28379 0012	MikroFIX™ for KiddieGAIT®/KiddieROCKER®	M	Black
28379 0013	MikroFIX™ for KiddieGAIT®/KiddieROCKER®	L/XL	Black
28752 0010	MikroFIX™ for 2.0 & 2 1/2 Allard AFOs	XS	
28752 0011	MikroFIX™ for 2.0 & 2 1/2 Allard AFOs	S	Black
28752 0013	MikroFIX™ for 2.0 & 2 1/2 Allard AFOs	M/L/XL	Black



## MIKROFIX SELF-ADHESIVE HOOK

MikroFIX™ rolls are the same ultra thin high strength hook material that is used to adhere the SoftKit to Allard AFO. Self-adhesive backing for easy application. Use to replace MikroFIX™ on Allard AFOs or any other orthosis fabrication when hook is required. Available in Black or White. Sold in packages of one 27 yard roll.




Item No.	Description	Size	Color
25040 0000	MikroFIX™ Hook Roll	2" x 27 yds	White
25041 0000	MikroFIX™ Hook Roll	2" x 27 yds	Black



## VERTICAL TIBIA PADS

Vertical tibia pads are designed to create an open channel for the tibial crest when extra pressure relief is necessary for optimum patient comfort. Vertical Tibia Pads fit all adult size Allard AFOs. Only available in Black. Sold in sets of 2 pieces.



Item No.	Description	Size	Color
28753 0011	Vertical tibia pads	XS/S	Black
28753 0013	Vertical tibia pads	M/L/XL	Black



## T-STRAP™ FOR ALLARD AFO

T-Strap™ pulls the talocrural and subtalar joints laterally to position them back underneath the tibia. Only available in Black. Choose Left or Right. Sold in packages of 1 piece.

*Quick and easy to apply:*

*Distal end attaches to bottom of ToeOFF® footplate – proximal “T” wraps around ankle and attaches to the Allard AFO lateral strut. Easily adjusted to customize to height of the patient’s medial malleolus. Plastazote inner core may be gently heated with dry heat to improve shape and fit. T-Strap is intended only for use in conjunction with an Allard AFO and a custom orthotic appropriate to meet the patient’s foot positioning requirements.*



Item No.	Description	Size	Color
28777 1000	T-Strap™, Left	One size	Black
28777 2000	T-Strap™, Right	One size	Black



## NON-SKID INTERFACE

Non-Skid Interface is designed to be used in between an Allard AFO and a plastic foot orthotic to prevent friction and abrasion of carbon composite footplates. Especially useful for pediatric patients that acquire sand inside of their shoes while playing. Only available in Black. Sold in packages of one meter rolls.

Item No.	Description	Size	Color
28411 0000	Non-Skid Interface Roll	5.5” x 39”	Black

## SOFT AFO STRAP ROLL

High-quality, double woven strapping material for any of your orthotic strapping needs. Slightly elasticized for increased flexibility. Only available in Black. Choose between two widths. Sold in packages of one 11 yard roll.

Item No.	Description	Size	Color
25060 0000	Soft AFO Strap Roll	1.5" x 11 yds	Black
25061 0000	Soft AFO Strap Roll	2" x 11 yds	Black



## FOOT ORTHOTIC SHELL

Proper foot positioning is crucial to optimize biomechanical function and maximize Allard AFO product durability. Foot orthotic shells are fabricated with EFM (Extra Firm Material) and Microcell Puff Interface, both are Ethylene Vinyl Acetate (EVA) foam. EVA is a shock absorber that is heat moldable, has excellent shape retention, and is abrasion resistant. Only available in Blue. Sold in packages of 2 pieces (left and right).

Item No.	Description	Size	Color
68372 0011	Pediatric	Small	Blue
68372 0012	Pediatric	Medium	Blue
68372 0013	Pediatric	Large	Blue
68372 0014	Pediatric	X-Large	Blue
68000 0010	Adult	X-Small	Blue
68000 0011	Adult	Small	Blue
68000 0012	Adult	Medium	Blue
68000 0013	Adult	Large	Blue
68000 0014	Adult	X-Large	Blue
68000 0015	Adult	XX-Large	Blue



## WEDGEKIT

Posting wedges are used for correcting or adjusting medial & lateral slopes. Using a WedgeKIT and documenting it is an excellent way to stay compliant when fitting an Allard AFO. Kit includes one lateral wedge and 3 heel wedges. All are pressure sensitive adhesive (PSA) Backed EVA foam, Shore A60 durometer. Sold in packages of 4 pieces.

Item No.	EVA Color	Dimension	Slope
28766 0000	Gray	3.875"	10 mm - 1 mm
	Chocolate	3.875"	7 mm - 1 mm
	Black	3.875"	5 mm - 1 mm
	Chocolate	7.75"	7 mm - 1 mm



## CENTRAL FABRICATION ACCESSORIES WARRANTY

Allard USA obligation under this warranty, is limited to the replacement of any defective part or parts for 90 days, starting at date of receipt.

It does not cover damage by tools, application of heat, or excessive bending/flexing of, or to any accessory or part thereof. Allard USA will not accept any product returns with biohazard contamination.

## RETURN POLICY

To help us serve you better and assure you receive appropriate credit, please call toll free 888-678-6548 to receive a Return Authorization (RA) number prior to shipping the return.

## DELIVERY TIMES

It is our goal to ship all orders received by 4:00pm, same day. If for any reason we are unable to meet that goal, you will be contacted.

## TERMS

There is no minimum order. There is a \$6.00 handling charge per order.

Payment terms are net 30 days. Overdue invoices are subject to a charge of 1.5% of the unpaid balance or .50% per month, whichever is greater.

## ORDERING INFORMATION

To place an order, please contact Customer Service:

Call toll free  
888-678-6548

Fax  
800-289-0809

Email  
orders@allardusa.com

## ADDRESS

Allard USA, Inc.  
300 Forge Way, Suite 3  
Rockaway, NJ 07866-2056

## HOURS OF OPERATION

Monday through Friday  
8:30AM - 6:00 PM EST

# Support for Better Life!

**allard** USA

allardusa.com

**ALLARD USA, INC.**  
300 Forge Way, Suite 3  
Rockaway, NJ 07866-2056

info@allardusa.com  
Toll Free 888-678-6548  
Fax 800-289-0809