

allard|USA

# X-LITE<sup>®</sup> PLUS

PLASTER MATERIAL for CASTING and SPLINTING

Say **goodbye to gloves**, masks and special ventilation. Choose an alternative solution **free from isocyanates** and **fiber glass**.

Protect our nature and **avoid** exposure to environmental and **health risks!**



## X-LITE® PLUS ROLLS & DISPENSERS

Our newest material for casting and splinting.  
 Rolls are available in 5 colors in 6 different widths.  
 Dispensers comes in Hot pink, Blue and Black in 15 cm widths.



Item No.	Color	Width	Length	Package
64701 0000	White	25 mm	1,8 m	5 rolls
64702 0000	White	50 mm	1,8 m	20 rolls
64703 0000	White	75 mm	1,8 m	20 rolls
64704 0000	White	100 mm	1,8 m	10 rolls
64706 0000	White	125 mm	1,8 m	10 rolls
64707 0000	White	150 mm	1,8 m	10 rolls
64711 0000	Red	25 mm	1,8 m	5 rolls
64712 0000	Red	50 mm	1,8 m	20 rolls
64713 0000	Red	75 mm	1,8 m	20 rolls
64714 0000	Red	100 mm	1,8 m	10 rolls
64716 0000	Red	125 mm	1,8 m	10 rolls
64717 0000	Red	150 mm	1,8 m	10 rolls
64721 0000	Blue	25 mm	1,8 m	5 rolls
64722 0000	Blue	50 mm	1,8 m	20 rolls
64723 0000	Blue	75 mm	1,8 m	20 rolls
64724 0000	Blue	100 mm	1,8 m	10 rolls
64726 0000	Blue	125 mm	1,8 m	10 rolls
64727 0000	Blue	150 mm	1,8 m	10 rolls
64731 0000	Black	25 mm	1,8 m	5 rolls
64732 0000	Black	50 mm	1,8 m	20 rolls
64733 0000	Black	75 mm	1,8 m	20 rolls
64734 0000	Black	100 mm	1,8 m	10 rolls
64736 0000	Black	125 mm	1,8 m	10 rolls
64737 0000	Black	150 mm	1,8 m	10 rolls
64741 0000	Hot Pink	25 mm	1,8 m	5 rolls
64742 0000	Hot Pink	50 mm	1,8 m	20 rolls
64743 0000	Hot Pink	75 mm	1,8 m	20 rolls
64744 0000	Hot Pink	100 mm	1,8 m	10 rolls
64746 0000	Hot Pink	125 mm	1,8 m	10 rolls
64747 0000	Hot Pink	150 mm	1,8 m	10 rolls
64750 0000	Mixed Color	150 mm	1,8 m	5 rolls
64791 0000	Blue	150 mm	10 m	1 dispenser
64792 0000	Black	150 mm	10 m	1 dispenser
64793 0000	Hot Pink	150 mm	10 m	1 dispenser

X-LITE® PLUS is suitable for all types of casts and splints. The material is made of a 100% cotton mesh, impregnated with a non-toxic thermoplastic resin. This provides you with a non-toxic and environmentally friendly material completely harmless to medical staff and patient. **EcoCast!**

- X-LITE® PLUS can be shaped easily and quickly on the patient without the use of gloves.
- X-LITE® PLUS is ultra-light and provides an airy and breathable cast.
- X-LITE® PLUS has a soft textile surface which provides additional comfort for the patient.



- X-LITE® PLUS provides a strong, durable cast with just 50% overlap.
- X-LITE® PLUS is re-usable and re-moldable any numbers of times.
- X-LITE® PLUS is available in 5 exciting colors – Black, White, Red, Blue and Hot Pink.



X-LITE® should be heated in a waterbath set at 70°C for 1-3 minutes. The material is fully cured after just 20 minutes.

X-LITE® PLUS can also be used for construction of removable splints and therapists can easily determine the stability of the splint by adding extra layers or individual reinforcements of the material. Hook, padding and other reinforcements attach easily and permanent on the X-LITE®.

X-LITE® PLUS can be re-heated and re-moulded several times for the same patient. Make small changes by spot heating a completed cast or splint instead of making a new one that will cost more, take more time and may cause further discomfort for your patient.

X-LITE® PLUS is available in rolls in 6 different widths and in a practical 15 cm width dispenser box.

## X-LITE® STOCKINETTE IN A PRACTICAL DISPENSER BOX!

We now offer a round-knitted Stockinette in a soft and durable quality.

The Stockinette is made from polyester and it shapes easily around the extremity, leaving no wrinkles against the skin. Available in a 25 meter roll dispenser box.



Item No.	Color	Width	Length	Package
64394 0000	White	30 mm	25 m	1 Dispenser
64395 0000	White	50 mm	25 m	1 Dispenser
64396 0000	White	70 mm	25 m	1 Dispenser
64397 0000	White	100 mm	25 m	1 Dispenser
64398 0000	White	150 mm	25 m	1 Dispenser



Learn more about how to work with our  
ECOCAST materials and find instructional  
Videos and Handbooks at [www.allardint.com](http://www.allardint.com).