

## Dorsiflexion Support for Foot Drop - Lightweight and Low Profile

NAVIGAIT™ A new alternative for patients with mild foot drop. The amount of dorsiflexion is easily adjusted with the elastic coil to accommodate the users specific needs.

- Supports ankle dorsiflexion during swing phase.
- Elastic coil adjustable and easily replaceable.
- Removable interface.
- The amount of dorsiflexion is adjusted through the position in the lacing and with the velcro attachment on the shin plate.
- Low profile design makes it invisible under slacks.
- Anterior shell can be trimmed with scissors if needed.

## NAVIGAIT™ & 4-FOOT Wrap

### **Intended use and Indications**

NAVIGAIT™ Intended to support the foot when the ability to actively dorsiflex is reduced due to conditions such as stroke, multiple sclerosis, post polio syndrome.

### **Contraindications**

NAVIGAIT™ should not be used by patients who present with spasticity or severe edema.


### **Limitations**

NAVIGAIT™ does not address malpositions of the foot or foot deformities. These should be cared for with shoes and/or orthotics.




The 4-FOOT Wrap makes it possible to wear NaviGAIT™ with sandals, shoes that do not have an anchor point, or barefoot. The wrap has a plastic insert for increased stabilization and non-slip silicon strips to keep it in place.

4-FOOT Wrap is available in 3 sizes. Measure circumference at metatarsal heads for correct size.



Item No.	Description	Size	Anterior section height
28410 0010	NAVIGAIT™	X-Small	23cm / 9"
28410 0011	NAVIGAIT™	Small/Medium	26.5cm / 10.5"
28410 0013	NAVIGAIT™	Large/XLarge	30cm / 11.75"
28412 0010	Additional Strap w/ Coil	X-Small	
28412 0011	Additional Strap w/ Coil	Small/Medium	
28412 0013	Additional Strap w/ Coil	Large/XLarge	
28413 0010	Additional Strong Strap w/ Coil	X-Small	
28413 0011	Additional Strong Strap w/ Coil	Small/Medium	
28413 0013	Additional Strong Strap w/ Coil	Large/XLarge	
28414 0010	Extra Strap Set	X-Small	
28414 0011	Extra Strap Set	Small/Medium	
28414 0013	Extra Strap Set	Large/XLarge	
			<b>Circumference</b>
28415 0011	4-FOOT Wrap	Small	(21-23 cm) 8.25" - 9"
28415 0012	4-FOOT Wrap	Medium	(23-26 cm) 9" - 10.25"
28415 0013	4-FOOT Wrap	Large	(26-29 cm) 10.25"-11.25"



Patent Pending

CampScandinavia\_Jan2019©