

GliderKIT™

Another Interface Innovation from Allard USA

Allard USA offers you yet another solution to meet the varying needs of your individual patients!



GliderKIT™

GliderKIT™ provides a gliding “sling” between the leg and the ToeOFF® anterior shell. It moves as the leg moves, making it ideal for sports activities such as running or bicycling. It may also be the interface of choice for individuals who continue to experience pressure on the tibia even though the ToeOFF is properly aligned, or for patients with small proximal tibia protuberances that press against the AFO anterior support.

- Cradles and moves with the leg to offer good feeling of support and comfort
- Especially while running or bicycling, reduces friction between leg and brace
- Breathable, thin, lightweight material



Moves with the leg, great for running or bicycling.

Creates a small “gap” between the leg and the brace to prevent pressure on the proximal tibia.



Select size that corresponds to ToeOFF Size

Select Left or Right based on desired D-Ring position		
Left Lateral Right Medial	Right Lateral Left Medial	Size
28355 1010	28355 2010	X-Small
28355 1011	28355 2011	Small
28355 1012	28355 2012	Medium
28355 1013	28355 2013	Large
28355 1014	28355 2014	X-Large