

More Allard AFOs with $\frac{1}{2}$ Heel Height!

ToeOFF® $2\frac{1}{2}$

BlueROCKER® $2\frac{1}{2}$

ToeOFF® $2\frac{1}{2}$ Addition





ToeOFF® 2½ and BlueROCKER® 2½

Same Patient and Practitioner friendly features you loved in ToeOFF® 2.0 and BlueROCKER® 2.0!

- Starter SoftKIT™ Included
- All MikroFIX Pre-Applied
- Shorter Wings - More Gently Contoured
- Alligator-type Tab System for Strap Application
- Choice of Wrap-Around or D-Ring Straps

½ Heel Height

- Easier to adapt to shoes with lower heel height
- Easier to customize to meet alignment and biomechanical needs
- Shifts anterior tibia plate alignment slightly forward
- Lower Curvature in the Fore foot - More space in the toe box

2½ Stability

- ToeOFF® 2½ - Same as ToeOFF® 2.0
- BlueROCKER® 2½ - Same as BlueROCKER® 2.0

	Ypsilon® Flow ½	ToeOFF® Flow 2½ KiddieFlow™	ToeOFF® 2½ ToeOFF® 2.0 ToeOFF® Short KiddieGAIT®	BlueROCKER® 2½ BlueROCKER® 2.0 BlueROCKER® Short KiddieROCKER®
Allows Range Of Motion	● ● ● ● ○	● ● ● ● ○	● ● ○ ○ ○ ○	● ○ ○ ○ ○ ○
M-L Stability	● ● ○ ○ ○ ○	● ● ● ● ○	● ● ● ● ● ○	● ● ● ● ● ●
A-P Stability	● ● ○ ○ ○ ○	● ● ● ● ○	● ● ● ● ● ○	● ● ● ● ● ●
Dorsiflexion Assist	● ● ○ ○ ○ ○	● ● ● ● ○	● ● ● ● ● ○	● ● ● ● ● ●
Spasticity Control	● ● ○ ○ ○ ○	● ● ● ● ○	● ● ● ● ● ○	● ● ● ● ● ●
Proximal Control	● ● ○ ○ ○ ○	● ● ● ● ○	● ● ● ● ● ○	● ● ● ● ● ●

*Approximate, varies slightly dependent on size.

Toe OFF® 2½ Addition

Three attractive surface options for your fashion conscious patients!

- Smooth, sleek surface allows pant legs to slide easily over brace
- Colored surface extends down the strut



CAMO

Stylish camouflage pattern.



BIRCH

Classic natural pattern blends with various wardrobe options.



BLACK

Sophisticated matte black provides a soft chic look and feel.

Size	ToeOFF® 2½		BlueROCKER® 2½		ToeOFF® 2½ Addition		
	w/ D-ring	w/o D-ring	w/D-ring	w/o D-ring	Black w/D-Ring	Camo w/D-ring	Birch w/D-Ring
XS-Left	28982 1010	28980 1010	28972 1010	28970 1010	28954 1010	28955 1010	28956 1010
XS-Right	28982 2010	28980 2010	28972 2010	28970 2010	28954 2010	28955 2010	28956 2010
S-Left	28982 1011	28980 1011	28972 1011	28970 1011	28954 1011	28955 1011	28956 1011
S-Right	28982 2011	28980 2011	28972 2011	28970 2011	28954 2011	28955 2011	28956 2011
M-Left	28982 1012	28980 1012	28972 1012	28970 1012	28954 1012	28955 1012	28956 1012
M-Right	28982 2012	28980 2012	28972 2012	28970 2012	28954 2012	28955 2012	28956 2012
L-Left	28982 1013	28980 1013	28972 1013	28970 1013	28954 1013	28955 1013	28956 1013
L-Right	28982 2013	28980 2013	28972 2013	28970 2013	28954 2013	28955 2013	28956 2013
XL-Left	28982 1014	28980 1014	28972 1014	28970 1014	28954 1014	28955 1014	28956 1014
XL-Right	28982 2014	28980 2014	28972 2014	28970 2014	28954 2014	28955 2014	28956 2014

Interface options for 2.0, 2½ and FLOW models of Allard AFOs

Because every TIBIA is different!



SoftKIT™ is the standard version. Soft foam in 6mm. **BLACK**



SoftKIT™ is the standard version. Soft foam in 6mm. **BEIGE**



ComfortKIT™ is the deluxe version of SoftKIT with 9mm memory foam. Interface of choice for patients with skin sensitivity and diabetic patients. **BLACK** with grey binding.



GliderKIT™ forms a sling that moves with the tibia. Recommended for high activity level patients, such as marathoners or cyclists. **BLACK**



CoverKIT™ uniquely combines a soft interface with a spandex sleeve to cover the brace. Great solution for fashion conscious patients. **BLACK**

Size table CoverKIT™

Size	Calf Circumference
X-Small	31 - 34 cm (12 ³ / ₁₆ - 13 ³ / ₁₆ in)
Small	34 - 37 cm (13 ³ / ₈ - 14 ⁹ / ₁₆ in)
Medium	37 - 41 cm (14 ¹ / ₁₆ - 16 ¹ / ₁₆ in)
Large	41 - 46 cm (16 ¹ / ₈ - 18 ¹ / ₈ in)
X-Large	46 - 51 cm (18 ¹ / ₈ - 20 ¹ / ₁₆ in)

Size table Interface Options

SoftKIT™ Black	SoftKIT™ Beige	ComfortKIT™	GliderKIT™ Left	GliderKIT™ Right	CoverKIT™ (size above)	Size of brace
28750 0010	28749 0010	28751 0010	28355 1010	28355 2010	28765 0010	X-Small
28750 0011	28749 0011	28751 0011	28355 1011	28355 2011	28765 0011	Small
28750 0012	28749 0012	28751 0012	28355 1012	28355 2012	28765 0012	Medium
28750 0013	28749 0013	28751 0013	28355 1013	28355 2013	28765 0013	Large
28750 0014	28749 0014	28751 0014	28355 1014	28355 2014	28765 0014	X-Large

allard

Allard USA Inc.

300 Forge Way, Suite 3, Rockaway NJ 07866-2056, USA | 888-678-6548 | info@allardusa.com | allardusa.com

Use only as directed. See IFU for full safety information, contraindications, and precautions.