



**More Allard AFOs  
with 1/2 Heel Height!**

ToeOFF<sup>®</sup> 2 1/2 addition



CAMO

BIRCH

BLACK

**NOW** More Options to Find the Best Functional Solution

# ToeOFF® 2 ½ BLUEROCKER® 2 ½



Same Patient and Practitioner friendly features you love in ToeOFF® 2.0 and BlueROCKER® 2.0!

- Starter SoftKIT™ Included
- All MikroFIX Pre-Applied
- Shorter Wings - More Gently Contoured
- Alligator-type Tab System for Strap Application
- Choice of Wrap-Around or D-Ring Straps

\*Approximate, varies slightly dependent on size.

## AND MORE!!!

- ½ Heel Height
  - Easier to adapt to shoes with lower heel height
  - Easier to customize to meet alignment and biomechanical needs
  - Shifts anterior tibia plate alignment slightly forward
- Lower Curvature in the Forefoot
  - More space in the toe box

**2 ½ STABILITY**  
 ToeOFF® 2 ½ - Same as ToeOFF® 2.0  
 BlueROCKER® 2 ½ - Same as BlueROCKER® 2.0

	Ypsilon® FLOW½	ToeOFF® FLOW2½ KiddieFLOW™	ToeOFF® 2½ ToeOFF® 2.0 ToeOFF® ToeOFF® Short KiddieGAIT™	BlueROCKER® 2½ BlueROCKER® 2.0 BlueROCKER® BlueROCKER® Short KiddieROCKER®
Allows Range Of Motion	●●●●○	●●●○○	●●○○○○	●○○○○
M-L Stability	●●○○○○	●●●○○	●●●●○	●●●●●
A-P Stability	●●○○○○	●●●○○	●●●●○	●●●●●
Dorsiflexion Assist	●●○○○○	●●●○○	●●●●○	●●●●●
Spasticity Control	●●○○○○	●●●○○	●●●●○	●●●●●
Proximal Control	●●○○○○	●●●○○	●●●●○	●●●●●

# INTERFACE OPTIONS!

For 2.0, 2 ½ and FLOW models of Allard AFOs

*Note to Patients: All padding options are available on Amazon.*

- SKINNY
- PROMINENT
- FLAT
- SENSITIVE
- WIDE
- IRREGULAR

every TIBIA is different!



SoftKIT™ is the standard version. Soft foam in 3mm. BLACK.



SoftKIT™ is the standard version. Soft foam in 3mm. BEIGE.



ComfortKIT™ is the deluxe version of SoftKIT with 5mm memory foam. Interface of choice for patients with skin sensitivity and diabetic patients. BLACK with gray binding.



CoverKIT™ uniquely combines a soft interface with a spandex sleeve to cover the brace. Great solution for fashion conscious patients. BLACK.

## Size table

SoftKIT™ Black	SoftKIT™ Beige	ComfortKIT™	CoverKIT™ (size below)	Size of brace
28750 0010	28749 0010	28751 0010	28765 0010	X-Small
28750 0011	28749 0011	28751 0011	28765 0011	Small
28750 0012	28749 0012	28751 0012	28765 0012	Medium
28750 0013	28749 0013	28751 0013	28765 0013	Large
28750 0014	28749 0014	28751 0014	28765 0014	X-Large

## Size table CoverKIT™

Size	Calf Circumference
X-Small	31 - 34 cm (12 3/16 – 13 3/8 in)
Small	34 - 37 cm (13 3/8 – 14 9/16 in)
Medium	37 - 41 cm (14 9/16 – 16 1/8 in)
Large	41 - 46 cm (16 1/8 – 18 1/8 in)
X-Large	46 - 51 cm (18 1/8 – 20 1/16 in)

# ToeOFF® 2 1/2 addition

Three attractive surface options for your fashion conscious patients:

- Smooth, sleek surface allows pant legs to slide easily over brace.
- Colored surface extends down the strut.



CAMO

**Stylish camouflage pattern.**



BIRCH

**Classic natural pattern blends with various wardrobe options.**



BLACK

**Sophisticated matte black provides a soft chic look and feel.**

AllardUSA\_AUG2023©

	ToeOFF® FLOW 2 1/2		ToeOFF® 2 1/2		BlueROCKER® 2 1/2		ToeOFF® 2 1/2 Addition		
	w/D-ring	w/o D-ring	w/D-ring	w/o D-ring	w/D-ring	w/o D-ring	Black w/D-Ring	Camo w/D-ring	Birch w/D-Ring
XS-Left	28962 1010	28960 1010	28982 1010	28980 1010	28972 1010	28970 1010	28954 1010	28955 1010	28956 1010
XS-Right	28962 2010	28960 2010	28982 2010	28980 2010	28972 2010	28970 2010	28954 2010	28955 2010	28956 2010
S-Left	28962 1011	28960 1011	28982 1011	28980 1011	28972 1011	28970 1011	28954 1011	28955 1011	28956 1011
S-Right	28962 2011	28960 2011	28982 2011	28980 2011	28972 2011	28970 2011	28954 2011	28955 2011	28956 2011
M-Left	28962 1012	28960 1012	28982 1012	28980 1012	28972 1012	28970 1012	28954 1012	28955 1012	28956 1012
M-Right	28962 2012	28960 2012	28982 2012	28980 2012	28972 2012	28970 2012	28954 2012	28955 2012	28956 2012
L-Left	28962 1013	28960 1013	28982 1013	28980 1013	28972 1013	28970 1013	28954 1013	28955 1013	28956 1013
L-Right	28962 2013	28960 2013	28982 2013	28980 2013	28972 2013	28970 2013	28954 2013	28955 2013	28956 2013
XL-Left	28962 1014	28960 1014	28982 1014	28980 1014	28972 1014	28970 1014	28954 1014	28955 1014	28956 1014
XL-Right	28962 2014	28960 2014	28982 2014	28980 2014	28972 2014	28970 2014	28954 2014	28955 2014	28956 2014

For more information on the 40 Reasons Allard became the world leader in carbon composite AFO technology, visit [AllardAFO.com](http://AllardAFO.com)!



ToeOFF® FLOW 2 1/2 info



BlueROCKER® 2 1/2 info

**Support for Better Life!**



[allardusa.com](http://allardusa.com)

**ALLARD USA, INC.**  
300 Forge Way, Suite 3  
Rockaway, NJ 07866-2056

[info@allardusa.com](mailto:info@allardusa.com)  
Toll Free 888-678-6548  
Fax 800-289-0809