

You Asked  
We Listened

More Allard AFOs with  
½ Heel Height!



ToeOFF<sup>®</sup> 2 ½ addition



NEW! CAMO

BIRCH

BLACK

**NOW** More Options to Find the Best Functional Solution

# ToeOFF® 2 ½ BLUEROCKER® 2 ½

Same Patient and Practitioner friendly features you love in ToeOFF® 2.0 and BlueROCKER® 2.0!

- Starter SoftKIT™ Included
- All MikroFIX Pre-Applied
- Shorter Wings - More Gently Contoured
- Alligator-type Tab System for Strap Application
- Choice of Wrap-Around or D-Ring Straps

## AND MORE!!!

- ½ Heel Height
  - Easier to adapt to shoes with lower heel height
  - Easier to customize to meet alignment and biomechanical needs
  - Shifts anterior tibia plate alignment slightly forward
- Lower Curvature in the Forefoot
  - More space in the toe box

15mm heel height

7mm heel height

\*Approximate, varies slightly dependent on size.

## 2 ½ STABILITY

ToeOFF® 2 ½ - Same as ToeOFF® 2.0  
BlueROCKER® 2 ½ - Same as BlueROCKER® 2.0


Call 888-678-6548 or e-mail [info@allardusa.com](mailto:info@allardusa.com) and ask for your FREE laminated Product Selection Guide!

Correct ASSESSMENT and CUSTOMIZATION is still important to give the patient optimal function from the Allard AFO.

Attention!

Please read Professional Instructions (included with each Allard AFO) regarding:

- Alignment
- Footplate trimming
- Correcting varus/valgus
- Shoe selection and Foot Orthotic Shell



2 1/2

Lower toe lift which results in more space in the toe box

You can also view our patient assessment and instructions on YouTube or receive CEUs by completing our online course!

ALLARD AFO FUNCTION SCALE

	Ypsilon® FLOW 1/2	ToeOFF® FLOW 2 1/2	ToeOFF® 2 1/2 ToeOFF® 2.0 ToeOFF® KiddieGAIT®	BlueROCKER® 2 1/2 BlueROCKER® 2.0 BlueROCKER® KiddieROCKER®
Allows Range Of Motion	● ● ● ● ○	● ● ● ○ ○	● ● ○ ○ ○	● ○ ○ ○ ○
M-L Stability	● ● ○ ○ ○	● ● ● ○ ○	● ● ● ● ○	● ● ● ● ●
A-P Stability	● ● ○ ○ ○	● ● ● ○ ○	● ● ● ● ○	● ● ● ● ●
Dorsiflexion Assist	● ● ○ ○ ○	● ● ● ○ ○	● ● ● ● ○	● ● ● ● ●
Spasticity Control	● ● ○ ○ ○	● ● ● ○ ○	● ● ● ● ○	● ● ● ● ●
Proximal Control	● ● ○ ○ ○	● ● ● ○ ○	● ● ● ● ○	● ● ● ● ●

# ToeOFF® 2 1/2 addition

Three attractive surface options for your fashion conscious patients:

- Smooth, sleek surface allows pant legs to slide easily over brace.
- Colored surface extends down the strut.



**NEW! CAMO**  
Stylish camouflage pattern.



**BIRCH**  
Classic natural pattern blends with various wardrobe options.



**BLACK**  
Sophisticated matte black provides a soft chic look and feel.

	ToeOFF® 2 1/2		BlueROCKER® 2 1/2		ToeOFF® 2 1/2 Addition		
	w/D-ring	w/o D-ring	w/D-ring	w/o D-ring	Black w/D-Ring	Camo w/D-ring	Birch w/D-Ring
XS-Left	28982 1010	28980 1010	28972 1010	28970 1010	28954 1010	28955 1010	28956 1010
XS-Right	28982 2010	28980 2010	28972 2010	28970 2010	28954 2010	28955 2010	28956 2010
S-Left	28982 1011	28980 1011	28972 1011	28970 1011	28954 1011	28955 1011	28956 1011
S-Right	28982 2011	28980 2011	28972 2011	28970 2011	28954 2011	28955 2011	28956 2011
M-Left	28982 1012	28980 1012	28972 1012	28970 1012	28954 1012	28955 1012	28956 1012
M-Right	28982 2012	28980 2012	28972 2012	28970 2012	28954 2012	28955 2012	28956 2012
L-Left	28982 1013	28980 1013	28972 1013	28970 1013	28954 1013	28955 1013	28956 1013
L-Right	28982 2013	28980 2013	28972 2013	28970 2013	28954 2013	28955 2013	28956 2013
XL-Left	28982 1014	28980 1014	28972 1014	28970 1014	28954 1014	28955 1014	28956 1014
XL-Right	28982 2014	28980 2014	28972 2014	28970 2014	28954 2014	28955 2014	28956 2014

For more information on the 40 Reasons Allard became the world leader in carbon composite AFO technology, visit [AllardAFO.com](http://AllardAFO.com)!



ToeOFF® 2 1/2 info



BlueROCKER® 2 1/2 info

Support for Better Life!