

F.O.S.

Foot Orthotic Shell



**Customize your ToeOFF® Family Product to meet specific patient needs**


- Optimize foot biomechanical function
- Enhance gait outcomes
- Maximize ToeOFF® product durability
- Assure patient safety

# F.O.S.

## Foot Orthotic Shell

### Description

A foot orthotic shell fabricated with EFM (Extra Firm Material) that is a very high density, heat moldable EVA. It is grindable and washable. It is the highest source of rigidity for orthoses, short of plastic. The Microcell Puff® Interface is Ethylene Vinyl Acetate (EVA) foam. EVA is a wonderful shock absorber that is heat moldable and has excellent shape retention and abrasion resistance.



Feature	Advantage	Benefit
EFM shell	<ul style="list-style-type: none"> <li>Heat moldable (275°)</li> <li>Easy to modify</li> <li>Easy to add posting</li> <li>Firm interface between patient and carbon footplate</li> </ul>	Optimize <ul style="list-style-type: none"> <li>Foot biomechanical function</li> <li>Gait outcomes</li> <li>ToeOFF product durability</li> <li>Patient safety</li> </ul>
Microcel Puff® interface	<ul style="list-style-type: none"> <li>Heat moldable</li> <li>Easy to modify</li> <li>Resistant to bottoming out</li> </ul>	Optimize <ul style="list-style-type: none"> <li>Patient comfort</li> <li>Orthotic shell durability</li> </ul>
	<ul style="list-style-type: none"> <li>Fit bilaterally</li> <li>Save cost of custom</li> </ul>	<ul style="list-style-type: none"> <li>No limb length discrepancy</li> <li>Added cash sale</li> </ul>

ORDERING INFORMATION	Adult	Description	Size	F.O.S Length	Ypsilon/ToeOFF/BlueRocker footplate length
	68000-0010	Foot Orthotic Shell,Adult XSmall, PR	XS	225mm	8 1/4" (210mm)
	68000-0011	Foot Orthotic Shell,Adult, Small, PR	S	240mm	9" (230mm)
	68000-0012	Foot Orthotic Shell,Adult Medium, PR	M	255mm	9 5/8" (245mm)
	68000-0013	Foot Orthotic Shell,Adult Large, PR	L	275mm	10 5/8" (270mm)
	68000-0014	Foot Orthotic Shell,Adult XLarge, PR	XL	290mm	11 1/4" (285mm)
	68000-0015	Foot Orthotic Shell,Adult XXLLarge, PR	XXL	305mm	N/A
	Pediatric				KiddieGAIT footplate length
	68372-0011	Foot Orthotic Shell, Pediatric Small, PR	S	150mm	6 1/4" (160mm)
	68372-0012	Foot Orthotic Shell, Pediatric Medium, PR	M	175mm	7 1/16" (180mm)
	68372-0013	Foot Orthotic Shell, Pediatric Large, PR	L	192mm	7 7/8" (200mm)
	68372-0014	Foot Orthotic Shell, Pediatric XLarge, PR	XL	210mm	N/A