

Give them the
stability they need
without locking up my ankles!

BENEFITS

- Functional heel-to-toe gait when used in conjunction with appropriate foot orthotic.
- Anterior support provides opportunity to enhance proprioceptive response.
- Lightweight—so important for these little feet.
- Can use less brace to do the same work as traditional plastic AFOs.
- Footplate gives spring assist - KiddieGAIT gives back energy.
- Anterior support adds stability and opportunity to manage knee extension/flexion moment.
- Softer and more dynamic heel rocker, which creates softer loading to reduce knee flexion moment.
- Therapy friendly - able to use SMO only for therapy, or with KiddieGAIT when stability and/or dynamic assist is required.
- Wearer compliance is increased - lightweight functional support allows more fluid and productive gait (less tiring).

KiddieGAIT™

KiddieROCKYER™

Available to fit infants 9 months of age to children 10 years of age!



KiddieGAIT offers a carbon composite framework for orthotist to build an orthosis to encourage functional heel-to-toe gait. While an SMO (molded inner boot) or UCB controls the foot, the KiddieGAIT floor reaction frame adds dynamic response, lower leg anterior support and stability, and/or opportunity for management of knee extension/flexion moment. This combination offers a unique opportunity for orthotist to customize function to meet specific patient needs.

Functional Heel-to-Toe Gait!



KiddieGAIT with UCB



KiddieROCKER with SMO



KiddieGAIT Baby Size with SureStep™

Patent no.:
ToeOFF, ToeOFF Short, ToeOFF Fantasy, ToeOFF NFR, BlueRocker, BlueRocker NFR, KiddieGAIT,
KiddieGAIT NFRAL736950, BE11005297, BE114626, DK1005297, DK114626, FI1005297,
FI114626, FR1005297, FR114626, IE1005297, IE114626, IT1005297, IT114626, CA2279225,
CNZL97181689, I, NL1005297, NL114626, NO313656, PL194247, CH1005297, CH114626,
ES1005297, ES114626, GB1005297, GB114626, SE1005297, SE114626, DE69709184-8-08
DE69732541-5-08, US5897515, ATE210417, ATE289187 Ypsilon, Ypsilon NFR, GE60208889-5-08,
GB1379201, IT1379201, US6887213

AllardUSA_MAR2015_A112