

WHY WILL  
THEY WANT  
TO WEAR?

- I love the “bounce in my step”.
- “So light” - not clunky!
- I can run and play with my friends and my sister.
- No pressure or blisters on my foot.
- I like that it looks so sporty!
- It is easy to put on - it stays in my shoe - I just slip my foot into my shoe.



**Trent Clayton**  
Allard “TeamUP”  
Junior Co-Captain.  
More info at:  
[getbackuptoday.com](http://getbackuptoday.com)


LET  
CHILDHOOD  
HAPPEN!

kiddieGAIT™  
kiddieROCKER™

YouTube

See kids in action!  
“Allard USA Pediatric Orthotics”:  
[www.youtube.com/allardusa](http://www.youtube.com/allardusa)

For more information  
contact your O&P provider!



AN ANKLE BRACE



THEY WILL  
WANT  
TO WEAR!



allard|USA  
Toll Free: 888-678-6548, Fax: 800-289-0809  
[www.allardusa.com](http://www.allardusa.com)

Find Us On Facebook & Twitter!



allard|USA  
[www.allardusa.com](http://www.allardusa.com)

Give them the **STABILITY** they need, without sacrificing... their **ABILITY** to play!



## MOBILITY BENEFITS

- Lightweight—so important for these little feet.
- Get the control that is needed with less brace.
- Anterior support adds stability and opportunity to manage knee hyperextension or crouch gait.
- Softer and more dynamic heel rocker, which creates softer loading to reduce knee flexion moment.
- Therapy friendly—able to use SMO only for therapy, or with KiddieGAIT and KiddieROCKER when stability and/or dynamic assist is required.
- It is no longer a matter of what they can or cannot do, but what they want to do.

## FUNCTIONAL HEEL-TO-TOE GAIT

KiddieGAIT and KiddieROCKER offer a carbon composite framework for your orthotist to build an orthosis to encourage functional heel-to-toe gait. While an SMO (molded inner boot) or UCB controls their foot, the KiddieGAIT and KiddieROCKER floor reaction frame adds dynamic response, lower leg anterior support and stability, and/or opportunity for management of knee extension/flexion moment. This combination offers a unique opportunity for your orthotist to customize function to meet specific needs of your child.



KiddieGAIT & KiddieROCKER address modern rehab and clinical techniques to achieve “as close to normal” gait as possible.

KiddieGAIT & KiddieROCKER should always be combined with an additional foot orthotic, designed to control the position of the foot.



Add a **fun** slip-on cover!

BABY SIZES

YOUTH SIZES



Note: Pattern & Colors are subject to change

KiddieGAIT and KiddieROCKER are available to fit infants 9 months of age to children 10 years of age. For older children, ask your orthotist about ToeOFF and BlueROCKER.

**allard** **USA**  
Support for Better Life!