

CUSTOM LENGTH UPRIGHTS

You asked, we listened!!!

For new adult patients and/or your current SWASH wearers that are growing up and current uprights just are not long enough, you may now custom order longer uprights.

- Add up to 6” distal length to any size
- Offered for 8mm diameter uprights only
- All SWASH Steady & Go except size 1 for Steady
- Old SWASH Low Profile

ORDERING INFORMATION

Item numbers:

288550000 SWASH STEADY & GO UPRIGHT SET 115° CUSTOM

288560000 SWASH STEADY & GO UPRIGHT SET 123° CUSTOM


Order specification

Please look at below illustration

1. Specify what size you need, based on Proximal length (A)
 - **MUST** be one of the six sizes indicated on the illustration
2. Specify desired additional (up to 6”) Distal Length (B) required

Allow 3-4 weeks for delivery.

Please send orders to orders@allardusa.com



Upright lengths

Size	Proximal Length A		Distal Length B		Total Length C	
	mm	in	mm	in	mm	in
1	63mm	2½in	175mm	7in	310mm	12¼in
2	80mm	3¼in	215mm	8½in	375mm	14¾in
3	80mm	3¼in	280mm	11in	440mm	17¼in
4	80mm	3¼in	320mm	12½in	480mm	19in
5	90mm	3½in	370mm	14½in	540mm	21¼in
6	100mm	4in	420mm	16½in	600mm	23½in

NOTE: If size selected in "Proximal Length A" doesn't provide clearance for the greater trochanter, go up one size.