

iFIT transtibial device

Fit and Aligned in a Single Session

Get those with limb loss
back on their feet quickly,
comfortably and affordably.



The iFIT transtibial device offers adjustability and customization that easily accommodates limb volume changes. Made of injection-molded, advanced polymer materials, it can be mass-produced cost effectively, benefiting those with insurance concerns and providing an affordable option for patients to get a secondary device.

- Immediate fit with no casts, scans or molds—customizable to fit patient needs
- Comfort & stability with customizable, padded liner
- High strength, lighter-weight socket
- Safety locking buckle system allows the user to adjust the socket size anytime
- Fit and aligned in a single visit with hand tools, ideal for any location
- Ideal as a preparatory prosthesis for the first year following amputation
- Perfect solution for patients with kidney or heart disease that experience daily volume changes
- Compatible with many commercially made feet and pylons

Adjust with iFIT!

To better serve you, we've partnered with Allard USA.



allardusa.com



“The best part about my iFIT leg is that I don’t have to worry about being in pain. It’s changed my life!”

ALLARD USA, INC.
300 Forge Way, Suite 3
Rockaway, NJ 07866-2056

info@allardusa.com
Toll Free 888-678-6548
Fax 800-289-0809



iFIT Prosthetics is a Veteran-Owned Small Business (VOSB)

Follow us: [f](#) [@](#) [You Tube](#)

iFIT transfemoral device

Adjust Anytime, Anywhere

Give those with
limb loss the mobility
they deserve.



iFIT transfemoral sockets have unique advantages for fitting persons with transfemoral limb loss including lower cost, adjustability, consistency of quality in manufacturing, and the immediacy of fit that speeds the rehabilitation process. A person can be fit in just under 2 hours.

- Buckle system easily accommodates highly compressible thigh tissue
- Light weight, strength tested, durable materials
- Customizable socket to maximize comfort for each individual patient
- Total contact giving a more stable gait and better range of motion
- Unique subischial design provides better sitting comfort
- Wide range of adjustments that can be made by the user
- Patients can adjust their fit any time without using socks
- A specially designed light weight offset offers rotation and various degrees of flexion of the socket to facilitate the alignment process

Adjust with iFIT!



“This prosthesis has been a godsend! There’s no words to describe how comfortable and easy it is.”

To better serve you, we’ve partnered with Allard USA.



The Right Fit, Right Now



Support for Better Life!

allardusa.com

ALLARD USA, INC.
300 Forge Way, Suite 3
Rockaway, NJ 07866-2056

info@allardusa.com
Toll Free 888-678-6548
Fax 800-289-0809



iFIT Prosthetics is a Veteran-Owned Small Business (VOSB)

Follow us: [f](#) [@](#) [You Tube](#)