



Help get patients back on their feet quickly

With iFIT, patients get faster, more efficient rehabilitation from acute care to outpatient gait training.

- Fit patients in just an hour with the iFIT transtibial adjustable prosthesis
- Control volume fluctuation throughout rehabilitation
- Easy adjustments can be made by therapist to maximize rehab outcomes
- Over the knee design aids in safe and stable gait pattern for those persons with diabetes or polyneuropathy
- Far more comfortable than traditional plaster casts with pylons and feet
- Perfect solution for patients with kidney, renal or heart disease that experience daily volume changes
- Therapists better control the rehabilitation process by optimizing the socket fit every day



iFIT PROSTHESIS FOR REHABILITATION AND AFTER

Distributed by

allardUSA

Support for Better Life!

Go to www.allardusa.com/iFIT-Prosthetics and watch a brief video on how the device works, or contact your Allard USA District Manager for more information.



iFIT Prosthetics is a Veteran-Owned Small Business (VOSB)



The Right Fit, Right Now

Follow us: [f](#) [@](#) [You Tube](#)