

COMBO Order Form

All Central Fabrication orders must be shipped to our Chattanooga location.

Date: _____

Contact: _____
P.O. #: _____
Phone: _____
Fax: _____
Email: _____

SHIP VIA: UPS FedEx Cust. P/U
SERVICE: Ground 3-Day 2-Day
 ND Saver Next Day ND Early AM

SHIP TO:	BILL TO:	<input type="radio"/> Same as Ship To
Acct #: _____	Acct #: _____	
Name: _____	Name: _____	
Address: _____	Address: _____	
Address 2: _____	Address 2: _____	
City, State, Zip Code: _____	City, State, Zip Code: _____	

PATIENT INFO

Patient ID: _____
Height: _____ Weight: _____

AFFECTED SIDE(S)

Left Knee Right Knee Both Knees


ALLARD AFO CHOICE

BlueROCKER[®] 2½ Recommended
 ToeOFF[®] 2½ May be considered Shoe Size: _____ (used as a guide to help determine Allard AFO size)

AFO PADDING OPTIONS

ComfortKIT[™] (Recommended) Premium 5mm memory foam interface for patients requiring added cushioning and a more intimate fit due to the shape of the tibia, skin condition, diabetes, and/or activity level.
 SoftKIT[™] Soft foam in 3mm.
 CoverKIT[™] Soft interface combined with a spandex sleeve to cover the AFO (Recommended for venous insufficiency)

AFO ALIGNMENT (TOE IN/OUT)


Left Foot		Right Foot
<input type="radio"/> Standard (5-7 degrees toe out) <input type="radio"/> Toe In: _____ (specify degrees) <input type="radio"/> Toe Out: _____ (specify degrees)		<input type="radio"/> Standard (5-7 degrees toe out) <input type="radio"/> Toe In: _____ (specify degrees) <input type="radio"/> Toe Out: _____ (specify degrees)

cont'd


PREPARE PATIENT FOR FITTING & LEG TRACING

Measurements are in: Inches CMs


1. Evaluate patient standing. Observe patients functional standing posture, particularly toe in/out position.




2. Patient should be supine. Be certain to duplicate standing and toe in/out posture observed in step one.




3a. Trace leg on tracing paper.



3b. Draw a horizontal line at level of Tibial Plateau.



3c. Take circumferential and M-L measurements, at knee center and 8" (20cm) above mid-patella.



SPECIAL INSTRUCTIONS:
