

VISSION™ Finger Splint POSI

A textile orthosis designed for clinical use in immobilization of the hand and fingers (II-V). The orthosis features an adjustable volar aluminum splint that enables precise positioning in the Position of Safe Immobilization (POSI) – a function that protects joints, reduces the risk of contractures, and promotes optimal healing. The adjustable distal strap can be positioned and customized to ensure an optimal fit. The wide strap helps reduce edema in the fingers and on the back of the hand.

- Stabilizes the wrist, palm, and all fingers
- Adjustable aluminum splint for customized fit
- Lightweight and soft, durable textile material for patient comfort
- Suitable for postoperative care, fractures, and rehabilitation

Recommended Range Of Application

Trauma, Post-Op, fractures of metacarpals and fingers. Dupuytren syndrome.

Measurements

Wrist Circumference

Item No.	Product Description	Size	Measurement
460191011	Vission Finger Splint POSI	Small	<16 cm
460191012	Vission Finger Splint POSI	Medium	16-19 cm
460191013	Vission Finger Splint POSI	Large	19-22 cm
460191014	Vission Finger Splint POSI	X-Large	>22 cm
460192011	Vission Finger Splint POSI	Small	<16 cm
460192012	Vission Finger Splint POSI	Medium	16-19 cm
460192013	Vission Finger Splint POSI	Large	19-22 cm
460192014	Vission Finger Splint POSI	X-Large	>22 cm

