

SELECTION Wrist Soft

Selection Wrist Soft is a standard wrist brace in a pull-on version with two closing tabs attached to lacing on the ulnar side of the wrist. It has an adjustable thumb strap, and the shape is low cut in the palm area to allow the hand to bend. This model is cut out around the thumb muscles to allow full movement of the thumb. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. It is intended to be used when the wrist needs support and off-loading in cases such as rheumatic disorders and soft tissue injuries.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Inflammation, work-related strain or when rheumatic and/or degenerative changes are present. Also for soft tissue injuries.

Contraindications

None known.

Measurements

Measure 1 = Wrist circumference

*Measure 2 = Circumference 8 cm (3.1") above the wrist. If this measurement is 18cm (7") or greater, choose the appropriate Plus size from the Article No. table below, otherwise choose the standard size based on Measure 1 (Wrist circumference).

This model is 18 cm (7") long on the volar side. Plus sizes have extra width in the forearm area.

SELECTION®



Item No.	Product Description	Size	L/R	Color	Measurement 1	Measurement 2
352041010	Selection Wrist Soft	X-Small	Left	Svart	< 11cm, 4.3in	< 11cm, 4.3in
352100000	-	-	-	-	22 - 168 mm	-
352041011	Selection Wrist Soft	Small	Left	Svart	11-14cm, 4.3in-5.5	11-14cm, 4.3in-5.5
352041115	Selection Wrist Soft	Small Plus	Left	Svart	11-14cm, 4.3in-5.5	18-21cm, 7.1in-8.27in
352041012	Selection Wrist Soft	Medium	Left	Svart	14-17cm, 5.5in-6.7in	14-17cm, 5.5in-6.7in
352041125	Selection Wrist Soft	Medium Plus	Left	Svart	14-17cm, 5.5in-6.7in	21-24cm, 8.27in-9.45in
352041013	Selection Wrist Soft	Large	Left	Svart	17-20cm, 6.7in-7.9in	17-20cm, 6.7in-7.9in
352041135	Selection Wrist Soft	Large Plus	Left	Svart	17-20cm, 6.7in-7.9in	24-27cm, 9.45in-10.63in
352041014	Selection Wrist Soft	X-Large	Left	Svart	20-23cm, 7.9in-9in	20-23cm, 7.9in-9in
352041145	Selection Wrist Soft	X-Large Plus	Left	Svart	20-23 cm, 7.9in-9in	27-30cm, 10.63in-11.8in
352041015	Selection Wrist Soft	XX-Large	Left	Svart	23-26cm, 9in-10.25in	23-26cm, 9in-10.25in
352041155	Selection Wrist Soft	XX-Large Plus	Left	Svart	23-26cm, 9in-10.25in	30-33cm, 11.1in-13in
352042010	Selection Wrist Soft	X-Small	Right	Svart	< 11cm, 4.3in	< 11cm, 4.3in
352042011	Selection Wrist Soft	Small	Right	Svart	11-14cm, 4.3in-5.5	11-14cm, 4.3in-5.5
352042115	Selection Wrist Soft	Small Plus	Right	Svart	11-14cm, 4.3in-5.5	18-21cm, 7.1in-8.27in
352042012	Selection Wrist Soft	Medium	Right	Svart	14-17cm, 5.5in-6.7in	14-17cm, 5.5in-6.7in
352042125	Selection Wrist Soft	Medium Plus	Right	Svart	14-17cm, 5.5in-6.7in	21-24cm, 8.27in-9.45in
352042013	Selection Wrist Soft	Large	Right	Svart	17-20cm, 6.7in-7.9in	17-20cm, 6.7in-7.9in
352042135	Selection Wrist Soft	Large Plus	Right	Svart	17-20cm, 6.7in-7.9in	24-27cm, 9.45in-10.63in
352042014	Selection Wrist Soft	X-Large	Right	Svart	20-23cm, 7.9in-9in	20-23cm, 7.9in-9in
352042145	Selection Wrist Soft	X-Large Plus	Right	Svart	20-23cm, 7.9in-9in	27-30cm, 10.63in-11.8in
352042015	Selection Wrist Soft	XX-Large	Right	Svart	23-26cm, 9in-10.25in	23-26cm, 9in-10.25in
352042155	Selection Wrist Soft	XX-Large Plus	Right	-	23-26cm, 9in-10.25in	30-33cm, 11.1in-13in