

# SELECTION Wrist Orthosis Rigid

SELECTION®



Selection Wrist Orthosis Rigid has the same features as the soft version but includes a thumb section for increased thumb support and stability. The thumb strap can be individually adjusted. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. This product is intended to be used when the first CMC-joint and wrist needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries or attritional wear.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

### Recommended Range Of Application

Inflammation, work related strain or when rheumatic and/or degenerative changes are present. Also, for soft tissue injuries and attritional wear.

### Contraindications

Non known.

### Measurements

Wrist circumference.

Item No.	Product Description	Size	Measurement	L/R	Color
352061010	Rigid	X-Small	< 11 cm	Left	Svart
352061011	Rigid	Small	11 - 14 cm	Left	Svart
352061012	Rigid	Medium	14 - 17 cm	Left	Svart
352061013	Rigid	Large	17 - 20 cm	Left	Svart
352061014	Rigid	X-Large	20 - 23 cm	Left	Svart
352061015	Rigid	XX-Large	23 - 26 cm	Left	Svart
352062010	Rigid	X-Small	< 11 cm	Right	Svart
352062011	Rigid	Small	11 - 14 cm	Right	Svart
352062012	Rigid	Medium	14 - 17 cm	Right	Svart
352062013	Rigid	Large	17 - 20 cm	Right	Svart
352062014	Rigid	X-Large	20 - 23 cm	Right	Svart
352062015	Rigid	XX-Large	23 - 26 cm	Right	Svart