

SELECTION Open Wrist

Selection Open Wrist is a brace that is simple to don and doff. It can be open totally flat. The integrated "wrap around" straps gives increased stability around the wrist. It has an adjustable thumb strap. The shape is cut out in the palm area which allows the hand to bend. This model is cut out around the thumb muscles to allow full movement of the thumb. It comes with a pre-shaped plastic spoon on the palmar side and a soft narrow plastic spoon on the dorsal side. It is intended to be used when the wrist needs off-loading in cases such as rheumatic disorders or after fracture healing.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Inflammation, work related strain or when rheumatic and/or degenerative changes are present.

Contraindications

Non known.

Measurements

Wrist circumference.

SELECTION®



Item No.	Product Description	Size	Measurement	L/R	Color
352071010	SELECTION Open Wrist	X-Small	< 11 cm	Left	Svart
352071011	SELECTION Open Wrist	Small	11 - 14 cm	Left	Svart
352071012	SELECTION Open Wrist	Medium	14 - 17 cm	Left	Svart
352071013	SELECTION Open Wrist	Large	17 - 20 cm	Left	Svart
352071014	SELECTION Open Wrist	X-Large	20 - 23 cm	Left	Svart
352072010	SELECTION Open Wrist	X-Small	< 11 cm	Right	Svart
352072011	SELECTION Open Wrist	Small	11 - 14 cm	Right	Svart
352072012	SELECTION Open Wrist	Medium	14 - 17 cm	Right	Svart
352072013	SELECTION Open Wrist	Large	17 - 20 cm	Right	Svart
352072014	SELECTION Open Wrist	X-Large	20 - 23 cm	Right	Svart