

# Selection Wrist Short

SELECTION Wrist Short is a flexible, elastic wrist orthosis with a textile reinforcement on the volar side. It is low cut in the thumb area to follow the hands natural bend and allow full movement of the thumb. Lacing closure makes it easy to don and doff with one hand. The lacing comes from the ulnar side so the risk of misalignment of the thumb is reduced. Thumb grip is adjustable. Can adjust the shape slightly by trimming with scissors. It is intended to be used when the wrist needs support and off-loading in cases such as rheumatic disorders, soft tissue injuries. Also, for attritional wear.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

## Recommended Range Of Application

Supports the wrist after injuries, operations, casting, hypermobility, rheumatic problems, chronic pain, osteoarthritis, overloading. Can also be used in preventive purpose during heavy or repetitive movements.

## Measurements

Wrist circumference.

SELECTION®



Item No.	Product Description	Size	Measurement	L/R	Color
352221010	Selection Wrist Short	X-Small	< 11 cm	Left	Svart
352221011	Selection Wrist Short	Small	11 - 14 cm	Left	Svart
352221012	Selection Wrist Short	Medium	14 - 17 cm	Left	Svart
352221013	Selection Wrist Short	Large	17 - 20 cm	Left	Svart
352221014	Selection Wrist Short	X-Large	20 - 23 cm	Left	Svart
352222010	Selection Wrist Short	X-Small	< 11 cm	Right	Svart
352222011	Selection Wrist Short	Small	11 - 14 cm	Right	Svart
352222012	Selection Wrist Short	Medium	14 - 17 cm	Right	Svart
352222013	Selection Wrist Short	Large	17 - 20 cm	Right	Svart
352222014	Selection Wrist Short	X-Large	20 - 23 cm	Right	Svart