

SELECTION Childrens Wrist

SELECTION®



Selection Childrens Wrist is a standard wrist brace in a, "glove" version with two closing tabs attached to lacing on the ulnar side of the wrist. It has an adjustable thumb strap and the shape is low cut in the palm area to allow the hand to bend. This model is also cut around the thumb muscles to allow full movement of the thumb. It comes with a simple aluminum stay.

It is intended to help guide the hand and wrist into a functional position. To reduce pain in combination with inflammatory and rheumatic conditions. To reduce spasticity in hand and wrist in neurological conditions.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Rheumatic conditions, neurological disorders where wrist position need mild correction or support, after casting for short period of use, and inflammatory conditions.

Contraindications

None known.

Measurements

Measure 1: Wrist circumference

Measure 2: Palm length

Item No.	Product Description	Size	L/R	Color	Length	Measurement 1	Measurement 2
353051009	SELECTION Childrens Wrist	XX-Small	Left	Svart	11 cm, 4.3in	8.5-11cm, 3.35in-4.3in	5-6cm, 2in-2.4in
353051010	SELECTION Childrens Wrist	X-Small	Left	Svart	13.5 cm, 5.3in	10-12.5cm, 3.95in-4.9in	6-7cm, 2.4in-2.75in
353051011	SELECTION Childrens Wrist	Small	Left	Svart	14.5cm, 5.7in	11.5-14cm, 4.5in-5.5in	7-8cm, 2.75in-3.15in
353051012	SELECTION Childrens Wrist	Medium	Left	Svart	17 cm, 6.7in	12-15.5cm, 4.7in-6.1in	8-9cm, 3.15in-3.54in
353051013	SELECTION Childrens Wrist	Large	Left	Svart	17.5 cm, 6.9in	14-16.5cm, 5.5in-6.5in	8-9cm, 3.15in-3.54in
353052009	SELECTION Childrens Wrist	XX-Small	Right	Svart	11 cm, 4.3in	8.5-11cm, 3.35in-4.3in	5-6cm, 2in-2.4in
353052010	SELECTION Childrens Wrist	X-Small	Right	Svart	13.5 cm, 5.3in	10-12.5cm, 3.95in-4.9in	6-7cm, 2.4in-2.75in
353052011	SELECTION Childrens Wrist	Small	Right	Svart	14.5cm, 5.7in	11.5-14cm, 4.5in-5.5in	7-8cm, 2.75in-3.15in
353052012	SELECTION Childrens Wrist	Medium	Right	Svart	17 cm, 6.7in	12-15.5cm, 4.7in-6.1in	8-9cm, 3.15in-3.54in
353052013	SELECTION Childrens Wrist	Large	Right	Svart	17.5 cm, 6.9in	14-16.5cm, 5.5in-6.5in	8-9cm, 3.15in-3.54in