

# NaviGAIT™ and 4-FOOT Wrap

NAVIGAIT™



NAVIGAIT™ is an alternative for users with mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the patient's specific needs.

Supports ankle dorsiflexion during swing phase. Elastic coil adjustable and easily replaceable. Contains removable interface. The amount of dorsiflexion is adjusted through the position in the lacing and with the velcro® attachment on the shin plate. Low profile design makes it invisible under slacks. Anterior shell can be trimmed with scissors if needed.

The optional 4-FOOT Wrap accessory makes it possible to wear NAVIGAIT™ with sandals, shoes that do not have an anchor point, or barefoot. The wrap has a plastic insert for increased stabilization and non-slip silicon strips to keep it in place.

## Recommended Range Of Application

Intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, or Spinal Cord Injury.

## Contraindications

NAVIGAIT™ should not be used by patients who present with spasticity or severe edema.

## Measurements

NAVIGAIT™: Height of plastic anterior section

NAVIGAIT™ 4-FOOT Wrap: Circumference of foot over the metatarsal heads

Item No.	Product Description	Size	Measurement	Height	Circumference
284100010	NAVIGAIT™	X-Small	-	9in / 23cm	-
284100011	NAVIGAIT™	Small/Medium	-	10.5in / 26.5cm	-
284100013	NAVIGAIT™	Large/X-Large	-	11.75in / 30cm	-
284150011	4-FOOT Wrap	Small	21-23 cm	-	8.25in - 9in / 21cm - 23cm
284150012	4-FOOT Wrap	Medium	23-26 cm	-	9in - 10.25in / 23cm - 26cm
284150013	4-FOOT Wrap	Large	26-29 cm	-	10.25in - 11.25in / 26cm - 29cm
284120011	Additional Strap w/ Coil	Small/Medium	-	-	-
284120013	Additional Strap w/ Coil	Large/X-Large	-	-	-
284130010	Additional Strong Strap w/ Coil	X-Small	-	-	-
284130011	Additional Strong Strap w/ Coil	Small/Medium	-	-	-
284140010	Extra Strap Set	X-Small	-	-	-
284140011	Extra Strap Set	Small/Medium	-	-	-
284140013	Extra Strap Set	Large/X-Large	-	-	-