

SoftKIT™

All foot drop 2.0 and 2½ Allard AFOs come with a Starter Interface Kit which includes one Starter SoftKIT™ and two Tibia Pads! Every ToeOFF® and BlueROCKER® should have a soft interface between the orthosis and the tibial crest. SoftKIT™ is a padding system that includes two narrow neoprene pads to create a pressure-free channel to help relieve any potential discomfort from bony prominences. An extra pad is included so one can be worn while the other is laundered.



Item No.	Product Description	Size	Color
283900010	SoftKIT™ for ToeOFF® & BlueROCKER®	X-Small	Black
287500010	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Small	Black
283910011	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Small	Black
283910012	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Medium	Black
283910013	SoftKIT™ for ToeOFF® & BlueROCKER® Short	Large	Black
283900011	SoftKIT™ for ToeOFF® & BlueROCKER®	Small	Black
287500011	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Small	Black
283900012	SoftKIT™ for ToeOFF® & BlueROCKER®	Medium	Black
287500012	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Medium	Black
283900013	SoftKIT™ for ToeOFF® & BlueROCKER®	Large/X-Large	Black
287500013	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Large	Black
287500014	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Large	Black
287490010	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Small	Beige
287490011	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Small	Beige
287490012	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Medium	Beige
287490013	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Large	Beige
287490014	SoftKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Large	Beige
283970008	SoftKIT™ for KiddieGAIT® Baby	BabySmall	Black
283970009	SoftKIT™ for KiddieGAIT® Baby	BabyMedium	Black
283970010	SoftKIT™ for KiddieGAIT® Baby	BabyLarge	Black
283780011	SoftKIT™ for KiddieGAIT® / KiddieROCKER® / KiddieFLOW™	Small	Black
283780012	SoftKIT™ for KiddieGAIT® / KiddieROCKER® / KiddieFLOW™	Medium	Black
283780013	SoftKIT™ for KiddieGAIT® / KiddieROCKER® / KiddieFLOW™	Large	Black
283780014	SoftKIT™ for KiddieGAIT® / KiddieROCKER® / KiddieFLOW™	X-Large	Black