


# MultiStatic Joint

The MultiStatic Joint enables fast, efficient, and effective static treatment of contractures due primarily to orthopedic injuries or diseases, utilizing HLBS (High Load Brief Stretch) principles. The angle of the joint is set by an eccentric axis. This continuously adjustable static extension or flexion joint is used in orthoses for both children and adults.

## Recommended Range Of Application

Improvement/rehabilitation of joint mobility and resting posture.



Item No.	Product Description	Size	Length	Diameter	Width	Height	Thickness	Weight
696890000	MultiStatic 9mm Static Joint	X-Small	100 mm	11 mm	9 mm	6.5 mm	2 mm	14 grams
696900000	MultiStatic 12mm Static Joint	Small	400 mm	19 mm	12 mm	8.5 mm	2 mm	83 grams
696910000	MultiStatic 16mm Static Joint	Medium	400 mm	24 mm	16 mm	9.5 mm	2 mm	122 grams
696920000	MultiStatic 19mm Static Joint	Large	400 mm	30 mm	19 mm	10 mm	2 mm	153 grams