

ToeOFF® Short

ToeOFF® Short is specifically developed to accommodate individuals with a shorter stature, slightly broader shaped tibia, and wider foot at the metatarsal area. Tibia padding should always be used for patient comfort. This model has a heel height of 7 mm in all sizes.

Recommended Range Of Application

Foot drop in combination with weak plantar flexors, and in some cases limited knee control. Mild to moderate spasticity and proximal deficits.

Contraindications

Should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits or severe spasticity.

ToeOFF® SHORT



| Item No. | Product Description | Size | L/R |
|-----------|-----------------------------------|--------|-------|
| 289901011 | ToeOFF® Short | Small | Left |
| 289901012 | ToeOFF® Short | Medium | Left |
| 289901013 | ToeOFF® Short | Large | Left |
| 289902011 | ToeOFF® Short | Small | Right |
| 289902012 | ToeOFF® Short | Medium | Right |
| 289902013 | ToeOFF® Short | Large | Right |
| 283910011 | ToeOFF® Short, SoftKIT™ | Small | - |
| 283910012 | ToeOFF® Short, SoftKIT™ | Medium | - |
| 283910013 | ToeOFF® Short, SoftKIT™ | Large | - |
| 287680011 | ToeOFF® Short, Wrap Around Straps | Small | - |
| 287680012 | ToeOFF® Short, Wrap Around Straps | Medium | - |
| 287680013 | ToeOFF® Short, Wrap Around Straps | Large | - |