

ToeOFF® Short

ToeOFF® Short is specifically developed to accommodate individuals with a shorter stature, slightly broader shaped tibia, and wider foot at the metatarsal area. Tibia padding should always be used for patient comfort. This model has a heel height of 7 mm in all sizes.

Recommended Range Of Application

Foot drop in combination with weak plantar flexors, and in some cases limited knee control. Mild to moderate spasticity and proximal deficits.

Contraindications

Should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits or severe spasticity.





Item No.	Product Description	Size
289901011	ToeOFF® Short	Small
289901012	ToeOFF® Short	Medium
289901013	ToeOFF® Short	Large
289902011	ToeOFF® Short	Small
289902012	ToeOFF® Short	Medium
289902013	ToeOFF® Short	Large
283910011	ToeOFF® Short, SoftKIT™	Small
283910012	ToeOFF® Short, SoftKIT™	Medium
283910013	ToeOFF® Short, SoftKIT™	Large
287680011	ToeOFF® Short, Wrap Around Straps	Small
287680012	ToeOFF® Short, Wrap Around Straps	Medium
287680013	ToeOFF® Short, Wrap Around Straps	Large