

F.O.S. Foot Orthotic Shells

Proper foot positioning is crucial to optimize biomechanical function and maximize Allard AFO product durability. Foot orthotic shells, FOS, are insoles fabricated with EFM (Extra Firm Material) and Microcell Puff® Interface, both are Ethylene Vinyl Acetate (EVA) foam. EVA is a wonderful shock absorber that is heat moldable and has excellent shape retention and abrasion resistance.



Item No.	Product Description	Size	Footplate Length	Length Insole
680000010	F.O.S. Adult, Pair	X-Small	8 1/4in	8 3/4in
680000011	F.O.S. Adult, Pair	Small	9in	9 1/2in
680000012	F.O.S. Adult, Pair	Medium	9 5/8in	10in
680000013	F.O.S. Adult, Pair	Large	10 5/8in	10 3/4in
680000014	F.O.S. Adult, Pair	X-Large	11 1/4in	11 1/2in
680000015	F.O.S. Adult, Pair	XX-Large	N/A	12in
683720011	F.O.S. Pediatric, Pair	Small	6 1/4in	5 3/4in
683720012	F.O.S. Pediatric, Pair	Medium	7 1/6in	6 3/4in
683720013	F.O.S. Pediatric, Pair	Large	7 7/8in	7 1/2in
683720014	F.O.S. Pediatric, Pair	X-Large	N/A	8 1/4in