

SoftSHELL™

SoftSHELL™ is a padding system that combines the SoftKIT™ interface with a lycra textile front to form a sleeve that slides over the anterior plate and a matching cover for the lower strut. SoftKIT™ is included.



Item No.	Product Description	Size	Color
283880008	SoftSHELL™ for KiddieGAIT® Baby	BabySmall	Blue
283880009	SoftSHELL™ for KiddieGAIT® Baby	BabyMedium	Blue
283880010	SoftSHELL™ for KiddieGAIT® Baby	BabyLarge	Blue
283890008	SoftSHELL™ for KiddieGAIT® Baby	BabySmall	Star Print
283890009	SoftSHELL™ for KiddieGAIT® Baby	BabyMedium	Star Print
283890010	SoftSHELL™ for KiddieGAIT® Baby	BabyLarge	Star Print
283590011	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Small	Skull Design
283590012	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Medium	Skull Design
283590013	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Large	Skull Design
283590014	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	X-Large	Skull Design
283580011	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Small	Pink Hearts
283580012	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Medium	Pink Hearts
283580013	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	Large	Pink Hearts
283580014	SoftSHELL™ for KiddieGAIT®, KiddieROCKER® & KiddieFLOW™	X-Large	Pink Hearts