

# Vertical Pads

Vertical pads are designed to create an open channel for the tibial crest when extra pressure relief is necessary for optimum patient comfort. They include pressure-sensitive backing for easy application. Intended as replacements for the originals supplied with each SoftKIT™.



Item No.	Product Description	Size	Measurement
287530011	Vertical Tibia Pads for 2.0 & 2½, Set of 2	XS/S	5/8in x 6 1/4in
287530013	Vertical Tibia Pads for 2.0 & 2½, Set of 2	M/L/XL	3/4in x 7 7/8in