

CoverKIT™

CoverKIT™ is another new smart interface innovation from Allard. This new interface is designed specifically for 2.0 and 2 ½ model AFOs. CoverKIT™ combines our SoftKIT™ interface with a thin textile front to form a sleeve that slides over the anterior shell of the AFO and gives it a smooth, sleek appearance. The kit includes two narrow neoprene pads if there is a need to create a pressure-free channel for tibia prominences. If needed, the original straps, supplied with the brace, can be worn together with the CoverKIT™.

Measurements

Mid-Calf Circumference



Item No.	Product Description	Size	Color	Calf circumference
287650010	CoverKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Small	Black	31–34 cm / 12 3/16–13 3/8 in
287650011	CoverKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Small	Black	34–37 cm / 13 3/8–14 9/16 in
287650012	CoverKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Medium	Black	37–41 cm / 14 9/16–16 1/18 in
287650013	CoverKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	Large	Black	41–46 cm / 16 1/8–18 1/8 in
287650014	CoverKIT™ for ToeOFF® & BlueROCKER® 2.0 & 2 ½	X-Large	Black	46–51 cm / 18 1/8–20 1/16 in