

# ComfortPAD™

The ComfortPAD™ provides a gliding "sling" between the leg and the Ypsilon® FLOW½ "wings" located on each side of the proximal "V" of the AFO. Easy to put on and take off, even with just one hand. Wide, thin band offers exceptional comfort with no "hot spots" of pressure on back of leg, while also offering good range of adjustability and maximum durability. Cradles and moves with the leg to offer good feeling of support and comfort. Reduces friction between leg and the AFO when walking, running, or bicycling.

## Recommended Range Of Application

Ideal for support activities such as running or bicycling. May also be the interface of choice for individuals who continue to experience pressure on the tibia even though their Ypsilon® FLOW½ is properly aligned, or for patients with small proximal tibia protuberances that press against the Ypsilon® FLOW½'s anterior support.



Item No.	Product Description	Size
288860000	Ypsilon® Flow ComfortPAD	One Size