

Sling with Abduction Pillow

Our Sling with Abduction Pillow transfers the weight of the forearm to the patient's back and to the unaffected shoulder while abducting the arm away from the body at about 10-15 degrees for shoulder rehabilitation.

Recommended Range Of Application

Rotator cuff repairs; Anterior repairs; Posterior dislocations; Capsular shifts; Global shoulder instability.

Measurements

Measurement: Length from elbow to wrist.



Item No.	Product Description	Size	Measurement
310230	Vission™ Sling with Abduction Pillow	Small	9in - 10in
310240	Vission™ Sling with Abduction Pillow	Medium	10in - 11in
310250	Vission™ Sling with Abduction Pillow	Large	11in - 12in
310260	Vission™ Sling with Abduction Pillow	X-Large	12in - 13in