

Vission™ Osteoarthritis Knee Orthosis

Our Vission™ Osteoarthritis Knee Orthosis offers an adjustable non-stretch strap that utilizes a 3-point counterforce system to apply force to the unaffected side of the knee, unloading the opposite affected compartment.

Measurements

Measurement: Knee circumference.

VISSION™



Item No.	Product Description	Size	Measurement	L/R
435941	Vission™ Osteoarthritis Knee Orthosis	Medium	14in - 15in	L Medial, R Lateral
435951	Vission™ Osteoarthritis Knee Orthosis	Large	15in - 17in	L Medial, R Lateral
435861	Vission™ Osteoarthritis Knee Orthosis	X-Large	17in - 19in	R Medial, L Lateral