

## CHECK™ Comfortable Hyper-Extension Control Knee

Soft upper and lower suspension sleeves and popliteal interface combine with adjustable straps to hold the frame in position as the straps are adjusted to provide a "stop" to control ability of the knee to hyper-extend.

## **Recommended Range Of Application**

Genu recurvatum, unstable knee caused by laxity of the ligaments, mild valgus instability, mild varus instability, neurologic function impairment (post C.V.A.).

## **Contraindications**

Not intended for sports activities

## Measurements

Posterior upright connection bar adjusts to fit knee widths from 3  $\frac{1}{2}$ " to 5  $\frac{1}{4}$ "

Item No.	Product Description	Measurement	Length	L/R
286520000	CHECK™ Stainless Steel	3½in - 5¼in	18in	-
286510000	Replacement Interface	N/A	N/A	Bilateral
286530000	D-Ring Set of 4	N/A	N/A	-



