

Vission™ LSO Short

Our Vission™ LSO Short combines rigid anterior and posterior panels by applying 360° of circumferential compression. By lifting the abdominal muscles, the brace supports and relieves stress on the lumbar spine.

Integrated rigid anterior and removable rigid posterior panels provide firm compressive stabilization. Anterior panels have finger loops for easy closure. The wide elastic cinch straps easily adjust to provide required compression to snug the posterior panel against spine. Rigid posterior panel is heat moldable. Ergonomic stays help to retain brace shape.

Recommended Range Of Application

Moderate to severe low back pain, lumbar sprains & strains, lumbar muscle weakness, lumbar instability.

Measurements

Measurement based on hip circumference.



Item No.	Product Description	Size	Measurement
201120	Vission™ LSO Short	X-Small	28.5in - 32.5in
201130	Vission™ LSO Short	Small	32.5in - 36.5in
201140	Vission™ LSO Short	Medium	36.5 - 40.5in
201150	Vission™ LSO Short	Large	40.5in - 44.5in
201160	Vission™ LSO Short	X-Large	44.5in - 48.5in
201170	Vission™ LSO Short	XX-Large	48.5in - 52.5in
201180	Vission™ LSO Short	3X-Large	52.5in - 60.5in
201190	Vission™ LSO Short	4X-Large	60.5in - 68.5in