

Vission™ LSO Tall

Our Vission™ LSO Tall combines rigid anterior and posterior panels by applying 360° of circumferential compression. By lifting the abdominal muscles, the brace supports and relieves stress on the lumbar spine.

Integrated rigid anterior and removable rigid posterior panels provide firm compressive stabilization. Anterior panels have finger loops for easy closure. The wide elastic cinch straps easily adjust to provide required compression to snug the posterior panel against spine. Rigid posterior panel is heat moldable. Ergonomic stays help to retain brace shape. Dual overlapping anterior panels for adjustable fit.

Measurements

Measurement based on hip circumference.



Item No.	Product Description	Size	Measurement
201220	Vission™ LSO Tall	X-Small	28.5in - 32.5in
201230	Vission™ LSO Tall	Small	32.5in - 36.5in
201240	Vission™ LSO Tall	Medium	36.5 - 40.5in
201270	Vission™ LSO Tall	XX-Large	48.5in - 52.5in