

Vission™ LSO+

The Vission™ LSO+ is now available to offer sagittal and coronal support in a lightweight, breathable back support. The back panel self molds to the anatomy forming an intimate fit. Combination of rigid front, side and back panels restricts lateral movement. All the benefits of our LSO Short with added sagittal and coronal support panels. Combine with our Extender to 'grow' orthosis 1-2 sizes larger.

Measurements

Measurement based on hip circumference.



Item No.	Product Description	Size	Measurement
203730	Vission™ LSO+	Small	32.5in - 36.5in
203740	Vission™ LSO+	Medium	36.5in - 40.5in
203760	Vission™ LSO+	X-Large	44.5in - 48.5