

# 3-Point Hyperextension Brace

Our Hyperextension Brace provides tri-planar stabilization in a lightweight easy to apply and adjust orthosis. Tri-planar control. Adjustable height. Spring loaded pelvic band reduces upward migration when sitting. Easy to use ratchet adjusts tension for back panel.

### Recommended Range Of Application

Moderate to severe low back pain, lumbar sprains & strains, lumbar muscle weakness, lumbar instability.

### Measurements

Measurement based on waist circumference.



Item No.	Product Description	Size	Measurement
200432	Vission™ 3-Point Hyperextension Brace	Small	24in - 29in
200442	Vission™ 3-Point Hyperextension Brace	Medium	29in - 35in
200452	Vission™ 3-Point Hyperextension Brace	Large	35in - 41in
200462	Vission™ 3-Point Hyperextension Brace	X-Large	41in - 45in