

Vission Thumb Orthosis Long with Stay

Soft and elastic wrist orthosis in a long model. The orthosis has a removable plastic stay that gives support to the thumb and surrounding area as well as light support for the wrist. This thumb orthosis is intended to restrict movements in CMC and MCP 1 joints and to position the thumb in a functional position. It is made of a durable and resilient material with a soft and comfortable lining. Open adjustable thumb grip.

Recommended Range Of Application

CMC and MCP 1 osteoartrosis or instability, Trigger Finger, Morbus De Quervain.



Item No.

460231009
460231010
460231011
460231012
460231013
460231014
460232009
460232010
460232011
460232012
460232013
460232014