

SELECTION Open Wrist Children

Selection Open Wrist Children is a brace that is simple to don and doff. It can be opened totally flat to make it easier to put on. The integrated "wrap around" straps give increased stability around the wrist. It has two thumb straps; one around the thumb and another on the dorsal side of the hand. This keeps the thumb in a correct position so that it does not drop towards the palm of the hand. The orthosis is intended to immobilize the CMC joint and wrist to reduce pain. It also guides the wrist into a neutral position to support a wrist after trauma and/or post operatively if needed.

Made with the airy and comfortable SELECTION material with all the functionality this incorporates.

Recommended Range Of Application

Rheumatoid Arthritis, Cerebral Palsy.

Measurements

Measure 1: Wrist circumference

Measure 2: Palm height

SELECTION®



Item No.	Product Description	Size	L/R	Color	Length	Measurement 1	Measurement 2
353061011	SELECTION Open Wrist Children	Small	Left	Svart	14.5 cm	11.5 - 14 cm	7 - 8 cm
353061012	SELECTION Open Wrist Children	Medium	Left	Svart	17 cm	12 - 15.5 cm	8 - 9 cm
353061013	SELECTION Open Wrist Children	Large	Left	Svart	17.5 cm	14 - 16.5 cm	8 - 9 cm
353062011	SELECTION Open Wrist Children	Small	Right	Svart	14.5 cm	11.5 - 14 cm	7 - 8 cm
353062012	SELECTION Open Wrist Children	Medium	Right	Svart	17 cm	12 - 15.5 cm	8 - 9 cm
353062013	SELECTION Open Wrist Children	Large	Right	Svart	17.5 cm	14 - 16.5 cm	8 - 9 cm