

# ToeOFF® FLOW 2½

For your active patients or others who require less resistance or more ROM. We applied Allard's 20+ years of AFO carbon composite experience and advances in technology to create a new proprietary formula to meet the functional and comfort needs of even more individuals.

The new FLOW technology offers a smoother transition (flow) throughout the gait cycle and has increased ROM in the sagittal plane. ToeOFF® FLOW 2½ has a lower heel height (approx. 7mm) which accommodates gentler contours of shoe insoles and offers more clearance for the forefoot in the shoe toe box. It is easier to customize to meet alignment and biomechanical needs of your patients.

## Recommended Range Of Application

Foot drop combined with 0 to moderate spasticity, in some cases even with impaired knee control and toe amputation. Limb proprioception deficit and mild proximal deficit. For further information on uses and contraindications see AFO Professional Instructions.

## Contraindications

ToeOFF® FLOW 2½ should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), or severe spasticity.

ToeOFF® FLOW 2½



Item No.	Product Description	Size	L/R
289601010	ToeOFF® FLOW 2½ w/o D-Ring	X-Small	Left
289601011	ToeOFF® FLOW 2½ w/o D-Ring	Small	Left
289601012	ToeOFF® FLOW 2½ w/o D-Ring	Medium	Left
289601013	ToeOFF® FLOW 2½ w/o D-Ring	Large	Left
289601014	ToeOFF® FLOW 2½ w/o D-Ring	X-Large	Left
289602010	ToeOFF® FLOW 2½ w/o D-Ring	X-Small	Right
289602011	ToeOFF® FLOW 2½ w/o D-Ring	Small	Right
289602012	ToeOFF® FLOW 2½ w/o D-Ring	Medium	Right
289602013	ToeOFF® FLOW 2½ w/o D-Ring	Large	Right
289602014	ToeOFF® FLOW 2½ w/o D-Ring	X-Large	Right
289621010	ToeOFF® FLOW 2½ with D-Ring	X-Small	Left
289621011	ToeOFF® FLOW 2½ with D-Ring	Small	Left
289621012	ToeOFF® FLOW 2½ with D-Ring	Medium	Left
289621013	ToeOFF® FLOW 2½ with D-Ring	Large	Left
289621014	ToeOFF® FLOW 2½ with D-Ring	X-Large	Left
289622010	ToeOFF® FLOW 2½ with D-Ring	X-Small	Right
289622011	ToeOFF® FLOW 2½ with D-Ring	Small	Right
289622012	ToeOFF® FLOW 2½ with D-Ring	Medium	Right
289622013	ToeOFF® FLOW 2½ with D-Ring	Large	Right
289622014	ToeOFF® FLOW 2½ with D-Ring	X-Large	Right