

ToeOFF® 2½ Addition

ToeOFF® 2½ Addition is a version of the ToeOFF® 2½ with a different color surface. The available options are Black, Birch or Camouflage. The surface gives the product a different look and design as well as a smoother surface to allow for clothing to slide more easy over the surface. This new manufacturing method also makes it possible to extend the colored surface all the way to the junction between the strut and the sole. Addition is only available with D-ring straps. Please note that the pattern in the Beige and Camouflage option may vary between orthoses, this is important to consider when ordering for bilateral wearers.

Recommended Range Of Application

ToeOFF® 2½ is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy, or Cerebral Palsy.

Contraindications

ToeOFF® 2½ should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

ToeOFF® 2½ addition



Item No.	Size	L/R	Product Description
289541010	X-Small	Left	ToeOFF® 2 ½ Addition, Black
289541011	Small	Left	ToeOFF® 2 ½ Addition, Black
289541012	Medium	Left	ToeOFF® 2 ½ Addition, Black
289541013	Large	Left	ToeOFF® 2 ½ Addition, Black
289541014	X-Large	Left	ToeOFF® 2 ½ Addition, Black
289542010	X-Small	Right	ToeOFF® 2 ½ Addition, Black
289542011	Small	Right	ToeOFF® 2 ½ Addition, Black
289542012	Medium	Right	ToeOFF® 2 ½ Addition, Black
289542013	Large	Right	ToeOFF® 2 ½ Addition, Black
289542014	X-Large	Right	ToeOFF® 2 ½ Addition, Black
289551010	X-Small	Left	ToeOFF® 2 ½ Addition, Camo
289551011	Small	Left	ToeOFF® 2 ½ Addition, Camo
289551012	Medium	Left	ToeOFF® 2 ½ Addition, Camo
289551013	Large	Left	ToeOFF® 2 ½ Addition, Camo
289551014	X-Large	Left	ToeOFF® 2 ½ Addition, Camo
289552010	X-Small	Right	ToeOFF® 2 ½ Addition, Camo
289552011	Small	Right	ToeOFF® 2 ½ Addition, Camo
289552012	Medium	Right	ToeOFF® 2 ½ Addition, Camo
289552013	Large	Right	ToeOFF® 2 ½ Addition, Camo
289552014	X-Large	Right	ToeOFF® 2 ½ Addition, Camo
289561010	X-Small	Left	ToeOFF® 2 ½ Addition, Birch
289561011	Small	Left	ToeOFF® 2 ½ Addition, Birch
289561012	Medium	Left	ToeOFF® 2 ½ Addition, Birch
289561013	Large	Left	ToeOFF® 2 ½ Addition, Birch
289561014	X-Large	Left	ToeOFF® 2 ½ Addition, Birch
289562010	X-Small	Right	ToeOFF® 2 ½ Addition, Birch
289562011	Small	Right	ToeOFF® 2 ½ Addition, Birch
289562012	Medium	Right	ToeOFF® 2 ½ Addition, Birch
289562013	Large	Right	ToeOFF® 2 ½ Addition, Birch
289562014	X-Large	Right	ToeOFF® 2 ½ Addition, Birch