

BlueROCKER® 2½

The Foot Drop Orthosis BlueROCKER® 2½ has the same dynamic functions and stability as a BlueROCKER® 2.0 with the exception that the 2½ has a lower heel height, 7 mm, to adapt to modern shoes. The toe lift is also slightly lower which results in more space in the toe box.

The BlueROCKER® 2½ has all MikroFIX™ pre-applied to simplify fitting of the padding. The strap solution is the same as on 2.0, made from a durable and strong material. The straps can be applied on the interior or external side of the AFO. It is always important to adapt any Allard AFO to the heel height of the shoe. Read more about this in the Allard AFO Professional Instructions.

Recommended Range Of Application

BlueROCKER® is designed to manage footdrop in conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, Spinal Cord Injuries, Traumatic Brain Injuries, Guillian-Barre Syndrome, Charcot-Marie-Tooth, Myelomeningocele, Neuropathy or Cerebral Palsy.

Contraindications

BlueROCKER® 2½ should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits (e.g. quadriceps spasticity, genu valgum or varum, genu recurvatum), severe spasticity.

BlueROCKER® 2½



Item No.	Product Description	Size	L/R
289701010	BlueROCKER® 2 ½ w/o D-Ring	X-Small	Left
289701011	BlueROCKER® 2 ½ w/o D-Ring	Small	Left
289701012	BlueROCKER® 2 ½ w/o D-Ring	Medium	Left
289701013	BlueROCKER® 2 ½ w/o D-Ring	Large	Left
289701014	BlueROCKER® 2 ½ w/o D-Ring	X-Large	Left
289702010	BlueROCKER® 2 ½ w/o D-Ring	X-Small	Right
289702011	BlueROCKER® 2 ½ w/o D-Ring	Small	Right
289702012	BlueROCKER® 2 ½ w/o D-Ring	Medium	Right
289702013	BlueROCKER® 2 ½ w/o D-Ring	Large	Right
289702014	BlueROCKER® 2 ½ w/o D-Ring	X-Large	Right
289721010	BlueROCKER® 2 ½ w/ D-Ring	X-Small	Left
289721011	BlueROCKER® 2 ½ w/ D-Ring	Small	Left
289721012	BlueROCKER® 2 ½ w/ D-Ring	Medium	Left
289721013	BlueROCKER® 2 ½ w/ D-Ring	Large	Left
289721014	BlueROCKER® 2 ½ w/ D-Ring	X-Large	Left
289722010	BlueROCKER® 2 ½ w/ D-Ring	X-Small	Right
289722011	BlueROCKER® 2 ½ w/ D-Ring	Small	Right
289722012	BlueROCKER® 2 ½ w/ D-Ring	Medium	Right
289722013	BlueROCKER® 2 ½ w/ D-Ring	Large	Right
289722014	BlueROCKER® 2 ½ w/ D-Ring	X-Large	Right